

2018 INTERNATIONAL PRODUCT CATALOG



drive

DeVilbiss
HEALTHCARE

TABLE OF CONTENTS

About Drive DeVilbiss Healthcare

Headquartered in Port Washington, New York, **Drive DeVilbiss Healthcare** manufactures a complete line of medical products, including mobility products, sleep and respiratory products, beds, bariatric products, wheelchairs, sleep surfaces and pressure prevention products, self-assist products, power operated wheelchairs, rehabilitation products, patient room equipment, personal care products and electrotherapy devices. Currently, the Company has corporate offices, manufacturing facilities and distribution facilities located throughout the United States, Canada, the United Kingdom, France, Germany, Romania, the Netherlands, China, Hong Kong, India and Australia. The Company markets its products to customers located throughout the United States, Canada, Mexico, Central and South America, Europe, the Middle East, Asia and Australia. For more information, www.drivemedical.com.

Growing our geographic reach through DDI

The new **Drive DeVilbiss International** activity unit was launched at the Medica show in Düsseldorf 2016. It is a key part of our strategy to grow global geographic reach. Staff from Drive Medical and DeVilbiss Healthcare form the new international team, increasing the company's manufacturing capabilities and its market presence around the world. The team is highly experienced in international commerce and with an exceptional product portfolio. DDI are the perfect partner for our international customers.

Drive DeVilbiss International provides sales, service and technical support from its offices in Mannheim, Germany and Tipton, UK. In addition, a number of product and sales specialists, who are located in key regions, support sales into Europe, Middle East, Africa and Asia. It also manages the rapidly growing business in India from an operations base set up in Mumbai, as well as an operational infrastructure in Sydney for Australia.

Healthcare requirements are rapidly developing and growing in many territories. The **Drive DeVilbiss International** team are reaching out to business to business distribution partners in these fast-paced markets to develop business together.

Mission Statement

As Drive DeVilbiss International our continuing mission is working with our partners to provide one of the widest ranges of hospital and home healthcare products available today.

The company is dedicated to its customers by consistently offering the best quality and the most innovative products across a broad spectrum of healthcare providers' needs, competitive pricing, education, world-class customer service, technical support and friendly people to work with.

Drive DeVilbiss Healthcare, enhancing the quality of life of the people that we touch.

Mobility	4
Walking Aids	18
Patient Room	34
Bath Safety	54
Kids	70
Respiratory	82



► **Mobility**

Wheelchairs

Transport Chairs

Shell Chairs

Comfort Chairs

Power Wheelchairs

Wheelchair Accessories

Blue Streak Wheelchair

This stylish wheelchair with blue powder-coated frame and reinforced cross brace comes standard with flip-back, padded desk-length armrests to help make transfers easier.

- Durable, lightweight, and easy-to-clean nylon upholstery
- Solid, durable, low-maintenance rubber tires mounted on composite wheels
- Push-to-lock wheel locks for added safe
- Calf straps for additional leg support



Easily Folded



ITEM #	SEAT WIDTH CM (IN)	ARMREST	PRODUCT FRONT RIGGING
BLS16FBD-SF	41 (16)	Flip Back Desk Arm	Swing-away Footrests
BLS18FBD-SF	46 (18)	Flip Back Desk Arm	Swing-away Footrests
BLS20FBD-SF	50 (20)	Flip Back Desk Arm	Swing-away Footrests

Specifications are available on page 15

Silver Sport 2 Wheelchair

This attractive silver vein steel-framed wheelchair is durable, stable and folds easily for storage and transport. Suitable for indoor and outdoor use, it features stylish 61 cm (24") Mag Wheels with solid, puncture-proof tires and resilient, lightweight plastic handrims.

- Comfortable, easy-to-clean padded armrests for support
- Easily accessible push-to-lock wheel locks on rear wheels to ensure safety
- Rear stepper tubes allow attendant to safely climb kerbs
- Adjustable tool-free swingaway footrests help facilitate easy transfers

ITEM #	SEAT WIDTH CM (IN)	ARMREST	PRODUCT FRONT RIGGING
SSP216FA-SF	41 (16)	Fixed Arm	Swing-away Footrests
SSP216FA-ELR	41 (16)	Fixed Arm	Swing-away, Elevating Leg rests
SSP218FA-SF	46 (18)	Fixed Arm	Swing-away Footrests
SSP218FA-ELR	46 (18)	Fixed Arm	Swing-away, Elevating Leg rests

Specifications are available on page 15

Stylish Mag Wheels Elevating Legrest



Stylish Mag Wheels Desk Armrests Easily Folded



Sentra EC Wheelchair

Designed to provide a comfortable, mobility experience for the bariatric user, this wheelchair with durable, carbon-steel frame and triple coated chrome provides an attractive, chip-proof finish for a polished look. The dual-axle position allows transition of seat height from standard to hemi level.

- Padded desk-style armrests are easily removable for transfer, transportation and storage
- Solid tires mounted on composite mag-style rear wheels provide durability and low maintenance
- 20 cm (8"), height-adjustable, solid front casters ensure durability
- Tool-free, height-adjustable footrests can be swung away to facilitate transfers

ITEM #	SEAT WIDTH CM (IN)	ARMREST	PRODUCT FRONT RIGGINGS
STD20ECDDA-SF	51 (20)	Detachable, Desk Arm	Swing-away, Adjustable Footrests
STD22ECDDA-SF	56 (22)	Detachable, Desk Arm	Swing-away, Adjustable Footrests
STD24ECDDA-SF	61 (24)	Detachable, Desk Arm	Swing-away, Adjustable Footrests
STD26DFA-SF	66 (26)	Detachable, Full Arm	Swing-away, Adjustable Footrests
STD28DFA-SF	71 (28)	Detachable, Full Arm	Swing-away, Adjustable Footrests
STD30DFA-SF	76 (30)	Detachable, Full Arm	Swing-away, Adjustable Footrests

Specifications are available on page 15

Sentra XXL Wheelchair

Designed with a 317.5 kg (700 lb.) weight capacity, this heavy-duty wheelchair offers reinforced steel gussets at all weight bearing points to provide additional strength. The riggings, sold separately, include extra-large, heavy-duty, aluminum foot plates for maximum support.

- Carbon steel frame with triple-coated chrome for an attractive, chip-proof, maintainable finish
- Heavy-duty, flat-free 20 cm (8") front casters— adjustable in two positions
- Push-to-lock wheel locks
- Extra-heavy-duty front forks
- Reinforced side frame, caster journals and dual cross bars



Extra-reinforced front wheel fork Reinforced dual cross brace Disc brakes optional




ITEM #	SEAT WIDTH CM (IN)	ARMREST	DISK BRAKE	PRODUCT FRONT RIGGINGS
STD28DDA-DB	71 (28)	Detachable Desk Arm	✓	Swing-away, Detachable Footrests
STD28DDA-SF	71 (28)	Detachable Desk Arm		Swing-away, Detachable Footrests
STD30DDA-DB	76 (30)	Detachable Desk Arm	✓	Swing-away, Detachable Footrests
STD30DDA-SF	76 (30)	Detachable Desk Arm		Swing-away, Detachable Footrests

Specifications are available on page 15



Crash Tested
in accordance with
ISO 7176-19

- Stepless adjustment of leg length 39–48 cm (15.4"–18.9")
- Seat height adjustment from 45–47.5 cm (17.7"–18.7") to in 3 steps via optional front wheels on the standard frame
- Removable and retractable leg supports
- Fold-up adjustable footrests adjustable in five steps from 0°–32°
- Standard quick-release mechanism for fast wheel mounting and dismounting
- Optional drum brake available



Rotec Standard Wheelchair

The Rotec Standard Wheelchair is available in eight seat widths, provides multiple adjustments and enables users to stay fit and active.

ITEM #	SEAT WIDTH CM (IN)	ARMREST	DRUM BRAKE	PRODUCT FRONT RIGGINGS
900100100	38 (15)	Adjustable Arm		Swing-away, Detachable Footrests
900100200	41 (16)	Adjustable Arm		Swing-away, Detachable Footrests
900100300	43 (17)	Adjustable Arm		Swing-away, Detachable Footrests
900100400	45 (18)	Adjustable Arm		Swing-away, Detachable Footrests
900100500	48 (19)	Adjustable Arm		Swing-away, Detachable Footrests
900100600	50 (20)	Adjustable Arm		Swing-away, Detachable Footrests
900100700	52 (20)	Adjustable Arm		Swing-away, Detachable Footrests
900100800	54 (21)	Adjustable Arm		Swing-away, Detachable Footrests
900150100	38 (15)	Adjustable Arm	✓	Swing-away, Detachable Footrests
900150200	41 (16)	Adjustable Arm	✓	Swing-away, Detachable Footrests
900150300	43 (17)	Adjustable Arm	✓	Swing-away, Detachable Footrests
900150400	45 (18)	Adjustable Arm	✓	Swing-away, Detachable Footrests
900150500	48 (19)	Adjustable Arm	✓	Swing-away, Detachable Footrests
900150600	50 (20)	Adjustable Arm	✓	Swing-away, Detachable Footrests
900150700	52 (20)	Adjustable Arm	✓	Swing-away, Detachable Footrests
900150800	54 (21)	Adjustable Arm	✓	Swing-away, Detachable Footrests

Specifications are available on page 15

Ecotec Standard Wheelchair

Available in five seat widths, the Ecotec Standard Wheelchair offers a cushioned backrest and seat cover for a more comfortable and customized fit for the user.

- Durable, corrosion-resistant steel frame
- Armrests can be rotated away for easy transfers
- Footrest can be folded up with the heel strap
- Six-position, height-adjustable leg support
- Optional drum brake available

ITEM #	SEAT WIDTH CM (IN)	ARMREST	DRUM BRAKE	PRODUCT FRONT RIGGINGS
900250100	39 (15)	Flip-up Arm		Swing-away Footrests
900250200	42 (17)	Flip-up Arm		Swing-away Footrests
900250300	45 (18)	Flip-up Arm		Swing-away Footrests
900250400	48 (19)	Flip-up Arm		Swing-away Footrests
900250500	51 (20)	Flip-up Arm		Swing-away Footrests
900200100	39 (15)	Flip-up Arm	✓	Swing-away Footrests
900200200	42 (17)	Flip-up Arm	✓	Swing-away Footrests
900200300	45 (18)	Flip-up Arm	✓	Swing-away Footrests
900200400	48 (19)	Flip-up Arm	✓	Swing-away Footrests
900200500	51 (20)	Flip-up Arm	✓	Swing-away Footrests

Specifications are available on page 15




Crash Tested
in accordance with
ISO 7176-19



Crash Tested
in accordance with
ISO 7176-19

- Handle height adjustable in four steps from 90–97.5 cm (35.5"–38.5")
- Back cushion adjustable in height and depth
- Three-position seat height adjustment in the standard model from 51–56 cm (20"–22") for maximum comfort
- Standard quick-release mechanism for fast wheel mounting and dismounting
- Stepless leg length adjustment 39–48 cm (15.4"–18.9")
- Heel straps provide additional safety



Freetec Lightweight Wheelchair

Available in eight seat widths, the Freetec Lightweight Wheelchair is highly customizable with three-position seat height, two-position armrest, and five-step footrest adjustments from 0°–32°. The handle height also adjusts in four steps and is available with drum brakes for added reliability.

ITEM #	SEAT WIDTH CM (IN)	ARMREST	DRUM BRAKE	PRODUCT FRONT RIGGINGS
900500100	38 (15)	Adjustable Arm		Adjustable Footrests
900500200	41 (16)	Adjustable Arm		Adjustable Footrests
900500300	43 (17)	Adjustable Arm		Adjustable Footrests
900500400	45 (18)	Adjustable Arm		Adjustable Footrests
900500500	48 (19)	Adjustable Arm		Adjustable Footrests
900500600	51 (20)	Adjustable Arm		Adjustable Footrests
900500700	52 (20)	Adjustable Arm		Adjustable Footrests
900500800	54 (21)	Adjustable Arm		Adjustable Footrests
900550100	38 (15)	Adjustable Arm	✓	Adjustable Footrests
900550200	41 (16)	Adjustable Arm	✓	Adjustable Footrests
900550300	43 (17)	Adjustable Arm	✓	Adjustable Footrests
900550400	45 (18)	Adjustable Arm	✓	Adjustable Footrests
900550500	48 (19)	Adjustable Arm	✓	Adjustable Footrests
900550600	51 (20)	Adjustable Arm	✓	Adjustable Footrests
900550700	52 (20)	Adjustable Arm	✓	Adjustable Footrests
900550800	54 (21)	Adjustable Arm	✓	Adjustable Footrests

Specifications are available on page 16

Wheelchairs
tested for transportation
of people in vehicles in
connection with all permit-
ted personal & wheelchair
restraint system



Crash Tested

XS2 Aluminum Wheelchair

Fitted with reflectors for outdoor visibility, the XS2 Aluminum Wheelchair comes standard with a seat cushion, tension-adjustable backrest and push-button length adjustable armrests. The simple folding mechanism, half-folding backrest, removable arms, quick-release footrests and self-propelling rear wheels allow for easy storage and transportation.

- Adjustable desk-style, flip-up armrests facilitate side transfers and give unhindered access to work surfaces
- Solid, durable, low-maintenance Polyurethane tires
- Rear stepper tubes with integrated anti-tip wheels for safe kerb climbing
- Push button lap belt for additional user stability

ITEM #	SEAT WIDTH CM (IN)	ARMREST	PRODUCT FRONT RIGGINGS	PRODUCT
XS2SP16SIL	40 (16)	Detachable Arm	Swing-away Footrests	Self Propel
XS2SP18SIL	45 (18)	Detachable Arm	Swing-away Footrests	Self Propel
XS2SP20SIL	50 (20)	Detachable Arm	Swing-away Footrests	Self Propel

Specifications are available on page 15



Three-position seat height adjustment



Ultra Lightweight Wheelchair

This Ultra Lightweight Wheelchair's aluminum frame is strong, stable and very lightweight for transporting and folding for storage. The half folding backrest, quick release footrests and rear wheels further facilitate storage and transportation.

- Attractive, stylish silver frame
- Comfortable, easy-to-clean, two-tone upholstery
- Soft touch fabric side panels improve comfort
- Desk style armrests allow easy access to work surfaces
- Convenient rear pocket for storing items

ITEM #	SEAT WIDTH (CM/IN)	ARMREST	PRODUCT FRONT RIGGINGS	PRODUCT
LAWC007A	62 (25)	Desk Arm	Quick-Release Footrests	Self Propel
LAWC011A	71 (28)	Desk Arm	Quick-Release Footrests	Self Propel

Specifications are available on page 16



Half Fold Back

Desk Style Armrests



Travelite Transport Chair

With a total weight of just 9,3 kg (20.5 lbs.), the lightweight Travelite Transport Chair is the ideal travel companion, whether traveling by car, plane or train. Stable, comfortable and able to fit through any doorway, this transport chair with 20 cm (8") transit wheels is ready to use and ensures maximum mobility.

- Soft seat and backrest provide comfort
- Convenient roller bag with wheels provides space for personal items
- Easily folds in just a few steps



Practical carrying bag

Carrying bag with wheels

Easily Folded

ITEM #	SEAT WIDTH CM (IN)
TC005	45 (18)

Specifications are available on page 16



Expedition Plus Transport Chair

Ideal for indoor or outdoor use, the Expedition Plus Transport Chair with 30 cm (12") transit wheels features a light aluminum frame for easy maneuverability. The convenient, easy-to-operate cable brakes allow for attendant to reduce speed or lock brakes for parking.

- Fold down backrest allows easy storage and transportation
- Full-length armrests with upholstered arm pads increase user comfort and are easy to clean
- Nylon 48 cm (19") seat and back offer comfortable support and are durable and easy to clean



Backrest can be folded up

Retractable leg supports

Foldable

ITEM #	SEAT WIDTH CM (IN)	COLOUR
EXP19LTBL	48 (19)	Blue

Specifications are available on page 16



Integra Shell Chair

Developed as an aid for moving residents around in nursing home environments, the Integra Shell Chair is useful for users who are seated for long periods at a time and at risk of developing pressure ulcers, providing additional support, comfort and positioning assistance.

- Tilt-in-space chair – gas strut operated
- Easy inclination due to the ergonomically designed push handle with hand lever
- Angle-adjustable and removable legrest
- Lap belt for additional safety
- Convenient storage pocket on the backrest
- Easy-to-manoeuvre, high-quality brake wheels (12.5 cm) which integrally rotate 360 degrees
- Good hoist access for easy moving and handling

ITEM #	SEAT WIDTH CM (IN)
INTEGRAM45CPI	45 (18)
INTEGRAM51CPI	51 (20)

Specifications are available on page 16



IDSoft Comfort Wheelchair

The IDSoft Comfort Wheelchair features an adjustable-height backrest and a headrest with multi-positional ear cushions. Available in three sturdy frame sizes, this wheelchair provides comfort and positioning to suit the user.

- Easily retractable legrests with adjustable calf cushions
- Removable, height- and depth-adjustable armrests
- Frame length adjustable according to the seat depth
- Width adjustment between legrests according to the seat depth

ITEM #	SEAT WIDTH CM (IN)	WHEELS	DRUM BRAKE
IDSFTSP16	39 (16)	Self Propel, 24" Wheels	
IDSFTSP18	44 (18)	Self Propel, 24" Wheels	
IDSFTSP20	49 (20)	Self Propel, 24" Wheels	
IDSFTTS16	39 (16)	Transit, 12" Wheels	
IDSFTTS18	44 (18)	Transit, 12" Wheels	
IDSFTTS20	49 (20)	Transit, 12" Wheels	
IDSFTMW16	39 (16)	Transit, 14.5" Midsize Wheels	
IDSFTMW18	44 (18)	Transit, 14.5" Midsize Wheels	
IDSFTMW20	49 (20)	Transit, 14.5" Midsize Wheels	
IDSFTSP16DB	39 (16)	Self Propel, 24" Wheels	✓
IDSFTSP18DB	44 (18)	Self Propel, 24" Wheels	✓
IDSFTSP20DB	49 (20)	Self Propel, 24" Wheels	✓



Powerstroll S Drive Wheelchair Pushing Aid

The heavy-duty Powerstroll S Drive Wheelchair Pushing Aid will turn most manual wheelchairs into caregiver-controlled powered wheelchairs, making pushing a wheel chair up hills or over troublesome gradients easier. Its unique, one-click docking station means the Powerstroll can be attached or removed easily in seconds.

- Allows speed of up to 4 mph
- Push-and-slide battery removal
- Hand control with variable speed dial, forward/reverse function
- 10-mile (16 km) travel on a full battery charge
- Fully programmable and comprehensive fault management system



ITEM #	CHARACTERISTICS
PWCPP009	Weight capacity (incl. chair) of 135 kg (298 lb)
PWCPP009HD	Bigger motor and weight capacity (incl. chair) of 180 kg (397 lb)





Wheelchair Poncho

Flat hem poncho with hood and zipper. Ideal for all wheelchair users.

- Nylon fabric with reflective strip for added safety
- Weather protection from unexpected showers
- Universal sizing with zip closure and drawstring hood

ITEM #	
RT-PONCHO1	



Wheelchair Backpack

The extremely spacious Backpack Bag is made from strong, durable and waterproof nylon, and is ideal for storing and carrying essentials when out and about.

- Convenient carry handle
- Reflective strips
- Waterproof nylon material provides protection
- Spacious, strong and durable

ITEM #	
RT-WCBAG1	



Mason Balance Air Cushion

The Mason Air Cushion provides an advanced level of pressure area care using a hand pump to adjust the cushion's air level to suit individual user comfort.

- Premier pressure management cushion
- Flexible interconnected, independent air cells
- Provides advanced level of pressure ulcer prevention
- Dynamic positioning and unmatched stability

ITEM #	SIZE CM (IN)
8047-1818-2	46x46 (18x18)



Vinyl Foam Cushion

- 5cm (2") thick foam
- Provides added comfort
- Velcro strips attached underneath cushion

ITEM #	SIZE CM (IN)
CU011	51x51 (20x20)
CU012	56x56 (22x22)

Item #	OVERALL DIMENSIONS CM (IN)				SEAT CM (IN)				WEIGHT KG (LBS)			UOM
	Overall Width Open	Overall Width Closed	Overall Length w/ Riggings	Overall Height	Seat Width	Seat Depth	Back Height	Seat to Floor Height	Product Weight	Carton Shipping Weight	Weight Capacity	
BLUE STREAK WHEELCHAIR												
BLS16FBD-SF	56 (22)	32 (12.5)	107 (42)	91 (36)	41 (16)	41 (16)	41 (16)	50 (19.5)	18 (40.25)	20 (45)	113.4 (250)	1/cs
BLS18FBD-SF	61 (24)	32 (12.5)	107 (42)	91 (36)	46 (18)	41 (16)	41 (16)	50 (19.5)	18.7 (41.25)	23.1 (51)	113.4 (250)	1/cs
BLS20FBD-SF	66 (26)	32 (12.5)	107 (42)	91 (36)	51 (20)	41 (16)	41 (16)	50 (19.5)	19 (42.25)	22.2 (49)	113.4 (250)	1/cs
SILVER SPORT 2 WHEELCHAIR												
SSP216FA-SF					41 (16)							1/cs
SSP216FA-ELR					41 (16)							1/cs
SSP218FA-SF					46 (18)							1/cs
SSP218FA-ELR					46 (18)							1/cs
SENTRA EC WHEELCHAIR												
STD20ECDDA-SF	71 (28)	32 (13)	107 (42)	91 (36)	51 (20)	46 (18)	41 (16)	43-48 (17-19)	22 (49)	29 (65)	204 (450)	1/cs
STD22ECDDA-SF	80 (31) 80	34 (13)	107 (42)	91 (36)	56 (22)	46 (18)	50 (20)	43-48 (17-19)	26 (57)		204 (450)	1/cs
STD24ECDDA-SF	85 (33) 80	37 (15)	107 (42)	91 (36)	61 (24)	46 (18)	50 (20)	43-48 (17-19)	26 (57)		204 (450)	1/cs
STD26DFA-SF	89 (35)	37 (15)	119 (47)	97 (38)	66 (26)	50 (20)	43-48 (17-19)	43-48 (17-19)	40 (88)		318 (701)	1/cs
STD28DFA-SF	94 (37)	37 (15)	119 (47)	97 (38)	71 (28)	50 (20)	43-48 (17-19)	43-48 (17-19)	40 (88)		318 (701)	1/cs
STD30DFA-SF	99 (39)	37 (15)	119 (47)	97 (38)	76 (30)	50 (20)	43-48 (17-19)	43-48 (17-19)	40 (88)		318 (701)	1/cs
SENTRA XXL WHEELCHAIR												
STD28DDA-DB	94 (37)	37 (14.5)	119 (47)	99 (39)	71 (28)	41 (16)	48-50	43-48 (17-19)	46,0 (101)	47 (103)	318 (700)	1/cs
STD28DDA-SF	94 (37)	37 (14.5)	119 (47)	99 (39)	71 (28)	41 (16)	48-50	43-48 (17-19)	46,0 (101)	47 (103)	318 (700)	1/cs
STD30DDA-DB	99 (39)	37 (14.5)	119 (47)	99 (39)	76 (30)	41 (16)	48-50	43-48 (17-19)	46,5 (103)	47 (104)	318 (700)	1/cs
STD30DDA-SF	94 (37)	37 (14.5)	119 (47)	99 (39)	76 (30)	41 (16)	48-50	43-48 (17-19)	46,5 (103)	47 (104)	318 (700)	1/cs
ROTEC STANDARD WHEELCHAIR												
900100100	57 (22)	28 (11)	106 (42)	96-103 (38-41)	38 (15)	43 (17)	51-56 (20-22)		20,2 (45)		125 (276)	1/cs
900100200	60 (24)	28 (11)	106 (42)	96-103 (38-41)	41 (17)	43 (17)	51-56 (20-22)		20,4 (45)		125 (276)	1/cs
900100300	62 (24)	28 (11)	106 (42)	96-103 (38-41)	43 (17)	43 (17)	51-56 (20-22)		20,5 (45)		125 (276)	1/cs
900100400	65 (26)	28 (11)	106 (42)	96-103 (38-41)	45 (18)	43 (17)	51-56 (20-22)		20,6 (45)		125 (276)	1/cs
900100500	68 (27)	28 (11)	106 (42)	96-103 (38-41)	48 (19)	43 (17)	51-56 (20-22)		20,7 (46)		125 (276)	1/cs
900100600	70 (28)	28 (11)	106 (42)	96-103 (38-41)	50 (20)	43 (17)	51-56 (20-22)		20,8 (46)		125 (276)	1/cs
900100700	72 (29)	28 (11)	106 (42)	96-103 (38-41)	52 (20)	43 (17)	51-56 (20-22)		20,8 (46)		125 (276)	1/cs
900100800	75 (30)	28 (11)	106 (42)	96-103 (38-41)	54 (21)	43 (17)	51-56 (20-22)		20,9 (46)		125 (276)	1/cs
900150100	57 (22)	28 (11)	106 (42)	96-103 (38-41)	38 (15)	43 (17)	51-56 (20-22)		20,2 (45)		125 (276)	1/cs
900150200	60 (24)	28 (11)	106 (42)	96-103 (38-41)	41 (17)	43 (17)	51-56 (20-22)		20,4 (45)		125 (276)	1/cs
900150300	62 (24)	28 (11)	106 (42)	96-103 (38-41)	43 (17)	43 (17)	51-56 (20-22)		20,5 (45)		125 (276)	1/cs
900150400	65 (26)	28 (11)	106 (42)	96-103 (38-41)	45 (18)	43 (17)	51-56 (20-22)		20,6 (45)		125 (276)	1/cs
900150500	68 (27)	28 (11)	106 (42)	96-103 (38-41)	48 (19)	43 (17)	51-56 (20-22)		20,7 (46)		125 (276)	1/cs
900150600	70 (28)	28 (11)	106 (42)	96-103 (38-41)	50 (20)	43 (17)	51-56 (20-22)		20,8 (46)		125 (276)	1/cs
900150700	72 (29)	28 (11)	106 (42)	96-103 (38-41)	52 (20)	43 (17)	51-56 (20-22)		20,9 (46)		125 (276)	1/cs
900150800	75 (30)	28 (11)	106 (42)	96-103 (38-41)	54 (21)	43 (17)	51-56 (20-22)		21 (46)		125 (276)	1/cs
XS2												
XS2SP16SIL	59 (23)	26 (10)	103 (40.5)	93 (37)	40 (16)	43 (17)	52 (21)		17 (37)		135 (298)	1/cs
XS2SP18SIL	63 (25)	26 (10)	103 (40.5)	93 (37)	45 (18)	43 (17)	52 (21)		17 (37)		135 (298)	1/cs
XS2SP20SIL	68 (27)	26 (10)	103 (40.5)	93 (37)	50 (20)	43 (17)	52 (21)		17 (37)		135 (298)	1/cs

Wheelchair Specifications

Item #	OVERALL DIMENSIONS CM (IN)				SEAT CM (IN)				WEIGHT KG (LBS)			UOM
	Overall Width Open	Overall Width Closed	Overall Length w/ Riggings	Overall Height	Seat Width	Seat Depth	Seat Height	Seat to Floor Height	Product Weight	Carton Shipping Weight	Weight Capacity	
FREETEC LIGHTWEIGHT WHEELCHAIR												
900500100	56 (22)		103 (41)	90-975 (35-38)	38 (15)	43-50 (17-20)	45-50 (18-20)		16 (32)		125 (276)	1/cs
900500200	59 (23)		103 (41)	90-975 (35-38)	41 (16)	43-50 (17-20)	45-50 (18-20)		16,3 (36)		125 (276)	1/cs
900500300	61 (24)		103 (41)	90-975 (35-38)	43 (17)	43-50 (17-20)	45-50 (18-20)		16,4 (36)		125 (276)	1/cs
900500400	63 (25)		103 (41)	90-975 (35-38)	45 (18)	43-50 (17-20)	45-50 (18-20)		16,5 (36)		125 (276)	1/cs
900500500	66 (26)		103 (41)	90-975 (35-38)	48 (19)	43-50 (17-20)	45-50 (18-20)		16,6 (37)		125 (276)	1/cs
900500600	68 (27)		103 (41)	90-975 (35-38)	50 (20)	43-50 (17-20)	45-50 (18-20)		16,8 (37)		125 (276)	1/cs
900500700	70 (28)		103 (41)	90-975 (35-38)	52 (20)	43-50 (17-20)	45-50 (18-20)		17 (37)		125 (276)	1/cs
900500800	73 (29)		103 (41)	90-975 (35-38)	54 (21)	43-50 (17-20)	45-50 (18-20)		17,2 (38)		125 (276)	1/cs
900550100	56 (22)		103 (41)	90-975 (35-38)	38 (15)	43-50 (17-20)	45-50 (18-20)		16 (32)		125 (276)	1/cs
900550200	59 (23)		103 (41)	90-975 (35-38)	41 (16)	43-50 (17-20)	45-50 (18-20)		16,3 (36)		125 (276)	1/cs
900550300	61 (24)		103 (41)	90-975 (35-38)	43 (17)	43-50 (17-20)	45-50 (18-20)		16,4 (36)		125 (276)	1/cs
900550400	63 (25)		103 (41)	90-975 (35-38)	45 (18)	43-50 (17-20)	45-50 (18-20)		16,5 (36)		125 (276)	1/cs
900550500	66 (26)		103 (41)	90-975 (35-38)	48 (19)	43-50 (17-20)	45-50 (18-20)		16,6 (37)		125 (276)	1/cs
900550600	68 (27)		103 (41)	90-975 (35-38)	50 (20)	43-50 (17-20)	45-50 (18-20)		16,8 (37)		125 (276)	1/cs
900550700	70 (28)		103 (41)	90-975 (35-38)	52 (20)	43-50 (17-20)	45-50 (18-20)		16,9 (37)		125 (276)	1/cs
900550800	73 (29)		103 (41)	90-975 (35-38)	54 (21)	43-50 (17-20)	45-50 (18-20)		16,9 (37)		125 (276)	1/cs
ECOTEC STANDARD WHEELCHAIR												
900250100	59,5 (23)		104 (41)	95 (37)	39 cm	42 (17)	47-51 (19-20)		17,2 (38)		130 (287)	1/cs
900250200	62 (24)		104 (41)	95 (37)	42 (17)	42 (17)	47-51 (19-20)		16,4 (36)		130 (287)	1/cs
900250300	65 (26)		104 (41)	95 (37)	45 (18)	42 (17)	47-51 (19-20)		17,5 (39)		130 (287)	1/cs
900250400	68 (27)		104 (41)	95 (37)	48 (19)	42 (17)	47-51 (19-20)		18,7 (41)		130 (287)	1/cs
900250500	71 (28)		104 (41)	95 (37)	51 (20)	42 (17)	47-51 (19-20)		19,8 (44)		130 (287)	1/cs
900200100	59,5 (23)		104 (41)	95 (37)	39 cm	42 (17)	47-51 (19-20)		17,2 (38)		130 (287)	1/cs
900200200	62,5 (24)		104 (41)	95 (37)	42 (17)	42 (17)	47-51 (19-20)		17,8 (39)		130 (287)	1/cs
900200300	65,5 (26)		104 (41)	95 (37)	45 (18)	42 (17)	47-51 (19-20)		18,9 (42)		130 (287)	1/cs
900200400	68,5 (27)		104 (41)	95 (37)	48 (19)	42 (17)	47-51 (19-20)		20 (44)		130 (287)	1/cs
900200500	71,5 (28)		104 (41)	95 (37)	51 (20)	42 (17)	47-51 (19-20)		21,2 (47)		130 (287)	1/cs

Transport Chairs Specifications

ULTRA LIGHTWEIGHT WHEELCHAIR												
LAWC007A	62 (24)	32 (13)	107 (42)	93 (37)	42 (17)	39 (15)	52 (20)		13 (29)		115 (254)	1/cs
LAWC011A	71 (28)	32 (13)	109 (43)	93 (37)	51 (20)	39 (15)	52 (20)		13 (29)		125 (276)	1/cs
TRAVELITE TRANSPORT CHAIR												
TC005	57 (22.5)	29 (11.5)	84 (33)	95 (37.5)	46 (18)	46 (18)	44 (17.5)	48 (19)	9 (20.5)	10,5 (23.2)	113,3 (250)	1/cs
EXPEDITION PLUS TRANSPORT CHAIR												
EXP19LTBL	56 (22)	25 (10)	93 (37)	96 (38)	48 (19)	42 (17)	49 (19)	48 (19)	12,0 (26.5)	10,8 (24)	176,9 (300)	1/cs

Shell Chairs Specifications

Item #	Type	OVERALL DIMENSIONS (CM/IN)			SEAT (CM/IN)			WEIGHT (KG/LBS)			UOM
		Overall Width	Overall Length	Overall Height	Seat Width	Seat Depth	Seat Height	Product Weight	Carton Shipping Weight	Weight Capacity	
INTEGRA SHELL CHAIR											
INTEGRAM45CPI	M	68 (26.8)-72 (28)	78 (31)	129 (51)-137 (54)	45 (18)	41-47 (16-17)		29 (64)		120 (256)	1/cs
INTEGRAM51CPI	M	68 (26.8)-72 (28)	78 (31)	129 (51)-137 (54)	51 (20)	41-47 (16-17)		29 (64)		120 (256)	1/cs





Walking Aids

Rollators

Walkers

Canes & Crutches



Migo 2G Standard Rollator

Featuring adjustable push handles that extend to 100 cm (39.3"), the Migo 2G Rollator is roughly 25 percent lighter than other rollators, making it ideal for travel and indoor and outdoor use.

- Tool-free assembly of the cane holder
- Delivered standard with basket, tray and cane holder
- One-hand folding mechanism
- Equipped with 50 mm (2") wide wheels, closed rims and adjustable hand brakes



ITEM# 723600000



Seat with integrated grip



Wide wheels with closed rim

ITEM #	PRODUCT	COLOR
723600000	Migo 2G	Silver

Specifications are available on page 27



ITEM# 721100100

Basket and tray



Dual-function brake



Cane holder



ITEM# 722100000

Safety release in red



Dual-function brake



Road Standard Rollator

Featuring height-adjustable, ergonomic hand grips, front and rear reflectors and dual-function brakes, the Road Rollator is ideal for travelers.

- Wheel diameter/width 200/37 mm (7.9"/1.5")
- Easy-to-install basket and tray
- Pull cord beneath seat for easy folding

ITEM# 721100100

Cristallo Standard Rollator

Equipped with adjustable hand brakes, multiple reflectors and an automatic safety lock to prevent folding, the Cristallo Rollator provides ideal visibility and safety for users.

- A one-hand folding mechanism for safe folding
- Features ergonomic height-adjustable hand grips
- Puncture-proof tire

ITEM# 722100000

ITEM #	PRODUCT	COLOR
721100100	Road	Silver
722100000	Christallo	Silver

Specifications are available on page 27



ITEM# R8RD-59

Easily folded



ITEM# R726BL

Dual-function Break



Backrest



Lightweight Aluminum Rollator

The lightweight aluminum Rollator features comfortable and ergonomic handgrips with loop cable brakes that are the ideal solution for those with low dexterity.

- Easy-to-use folding mechanism
- 18 cm (7") diameter solid tires
- Curved flip-up backrest for comfort and easy storage

ITEM# R8RD-59 / R8BL-59

Aluminum Rollator

Equipped with 15.2 cm (6") non-marring casters and serrated-edge brakes, the Aluminum Rollator provides a smooth and safe ride whether indoors or outdoors.

- Seamless padded seat flips up for additional storage
- Removeable back rest that can fold up or down
- Easy-to-use deluxe loop locks
- Zippered pouch under seat for added privacy and security of personal items

ITEM# R726BL / R726RD


ITEM #	PRODUCT	COLOR
R8RD-59	Leightweight Aluminum	Red
R8BL-59	Leightweight Aluminum	Blue
R726BL	Aluminum Rollator	Blue
R726RD	Aluminum Rollator	Red

Specifications are available on page 27



ITEM# 745R

Wheel easily attached



Does not move when folded



ITEM# 770100100

Large seat



With backrest



Basket



Diamond Deluxe Rollator

Capable of converting from a rollator to a transport chair, the lightweight Diamond Deluxe Rollator adjusts to accomodate a range of individual needs. This rollator also features retractable and removable footrests.

- Wheel diameter/width 200/30 mm (7.9"/1.5")
- Height-adjustable hand grips
- One-step pull cord for easy folding

ITEM# 745R

GoLite 200 XXL Rollator

Featuring puncture-proof tires, ergonomic hand grips that adjust to six different positions, a cushioned, removable back support bar and a 200 kg (441 lb.) weight limit, the GoLite 200 XXL Rollator is sturdy and safe.

- Latched seat prevents folding
- Equipped with a padded, extra large seat
- Wheel diameter/width 195/50 mm (7.8"/2")
- Easy-to-set brake

ITEM# 770100100

ITEM #	PRODUCT	COLOR
745R	Diamond Deluxe	Red
745SI	Diamond Deluxe	Silver
770100100	GoLite 200 XXL	Red

Specifications are available on page XX



ITEM# RTL10266

Nitro Rollator


Equipped with encased brake wires and cross bar frames, the Nitro Rollator offers optimal stability and safety. This rollator also features extra large 25 cm (10") casters for idea steering and maneuverability.

- A one-hand folding mechanism for easy travel and storage
- Equipped with cane holder

ITEM# RTL10266 / RTL10266-BLK



ITEM# RTL10266CF



Nitro Carbon Rollator

Weighing 5.5 kg (12 lbs.), the Nitro Carbon Fiber Rollator features a high-quality, leather seat, cushioned back rest and ergonomic, height-adjustable hand grips for added comfort.

- Carbon fiber frame
- Extra-large 25 cm (10") casters
- Easy-to-fold design ideal for storage and travel
- Standard passive lighting

ITEM# RTL10266CF

ITEM #	PRODUCT	COLOR
RTL10266	Nitro	Red
RTL10266-BLK	Nitro	Black
RTL10266CF	Nitro Carbon	Black

Specifications are available on page 27

Roomba Rollator

Capable of navigating narrow hallways and tight corners, the Roomba Rollator is specially designed for indoor use. This rollator is equipped with ergonomic handles and an end-to-end brake lever that can be engaged with one hand.

- High-quality craftsmanship and materials
- Includes basket and tray

ITEM# 730200100


Ultralight Tri Walker Rollator

Weighing 4.5 kg (10 lbs.), this ultra-lightweight rollator features adjustable handles, 17.8 cm (7") front swivel wheels and an easy, one-hand folding capability for ideal mobility. This rollator is equipped with looped cable brakes that activate easily and lock for safety.


- One-handed folding mechanism for easy travel and storage
- Includes carry pouch
- Wire basket and tray are available as optional extras



ITEM# 730200100



ITEM# TW008B-23



ITEM #	PRODUCT	COLOR
730200100	Roomba	White
TW008B-23	Ultralight Tri Walker	Blue

Specifications are available on page 27



Item #	OVERALL DIMENSIONS CM (IN)				SEAT CM (IN)				WEIGHT KG (LBS)			UOM	Color
	Overall Width Open	Overall Width Closed	Overall Length	Overall Height	Seat Width	Seat Depth	Seat Height	Seat to Floor Height	Product Weight	Carton Shipping Weight	Weight Capacity		
MIGO 2G STANDARD ROLLATOR													
723600000	58 (23)	23 (9)	61 (24)	80-100 (31-39)	36 (14)	23,5 (9)	60 (24)		9.0 (19.8)		130 (287)	2/cs	Silver
ROAD STANDARD ROLLATOR													
721100100	59 (23)	25 (10)	61 (24)	91,5 (36)	39 (15)	23 (9)	60 (24)		9.6 (21.2)		120 (265)	1/cs	Silver
CRISTALLO STANDARD ROLLATOR													
722100000	61 (24)	26 (10)	67 (26)	10 (39)	38 (15)	16 (6)	49 (19)		10.3 (22.7)		130 (287)	1/cs	Silver
LIGHTWEIGHT ALUMINUM ROLLATOR													
R8RD-59	64 (25)	20 (8)	62 (24)	83-98 (33 - 38)	36 (14)	16 (6)	57 (24)		7.0 (15.4)		135 (298)	1/cs	Red
R8BL-59	64 (25)	20 (8)	62 (24)	83-98 (33 - 38)	36 (14)	16 (6)	57 (24)		7.0 (15.4)		135 (298)	1/cs	Blue
ALUMINUM ROLLATOR													
R726BL	62 (24)		70 (28)		36 (14)	31 (12)	51 (20)		7.2 (15.9)		136 (300)	1/cs	Blue
R726RD	62 (24)		70 (28)		36 (14)	31 (12)	51 (20)		7.2 (15.9)		136 (300)	1/cs	Red
DIAMOND DELUXE ROLLATOR													
745R	64 (25)	28 (11)	104,2 (41)	77-94,5 (30-37)	40 (16)	19 (7)	54 (21)		10.3 (22.7)		136 (300)	1/cs	Red
745SI	64 (25)	28 (11)	104,2 (41)	77-94,5 (30-37)	40 (16)	19 (7)	54 (21)		10.3 (22.7)		136 (300)	1/cs	Silver
GOLITE 200 XXL ROLLATOR													
770100100	64 (25)	28 (11)	104,2 (41)	77-94,5 (30-37)	40 (16)	19 (7)	54 (21)		10.3 (22.7)		136 (300)	1/cs	Red
NITRO ROLLATOR													
RTL10266	60,5 (24)	25 (10)	70,5 (28)	82-95 (32-37)	44,5 (17)	22 (9)	49 (19)		7.0 (15.4)		130 (287)	1/cs	Red
RTL10266-BLK	60,5 (24)	25 (10)	70,5 (28)	82-95 (32-37)	44,5 (17)	22 (9)	49 (19)		7.0 (15.4)		130 (287)	1/cs	Black
NITRO CARBON ROLLATOR													
RTL10266CF	59 (23.2)	24 (9.4)	70,5 (27.8)	84,5 (33.2)-97 (38.2)	44 (17.3)	22 (8.7)	53 (20.8)		5.8 (12.8)		150 (330.7)	1/cs	Black
ULTRALIGHT TRI WALKER ROLLATOR													
TW008B-23	61 (24)	27 (11)	56 (22)	78-88 (31-35)	44 (17.3)	22 (8.7)	53 (20.8)		5 (11)		5,0 (11.0)	1/cs	Blue
ROOMBA ROLLATOR													
730200100	57 (22)	29 (11)	55 (22)	81,4-94,1 (32-37)					5.6 (12.3)		100 (220)	1/cs	White



Front Wheels set to swivel position



Front Wheels set to fixed position



ITEM# 10230

Clever-Lite Walker - 5" Casters

With front wheels capable of changing from a swivel to fixed position at any time, the Clever-Lite Walker allows the user to adjust the walker to meet their individual needs.

- Sturdy 2.5 cm (1") diameter anodized aluminum construction
- Includes flip-up seat with built-in carry handle
- Choice of 12.7 cm (5")
- Sophisticated loop locks and patented dual-lever side paddle release

ITEM# 10230

Folding Walking Frame

The lightweight Folding Walking Frame is ideal for users with varying degrees of hand strength and dexterity. This walker is equipped with a unique "U" frame design that allows the user to activate the individually operating sides by pressing an easy-push button with their fingers, palms or side of hands.

- Height-adjustable legs
- Available with or without wheels
- Perfect for sit-to-stand assistance
- Easy to fold for storage and travel

ITEM# WA008 / WA009



ITEM# WA003

Front Wheels set to fixed position



ITEM# 700100100

Folding



Adjustable legs



Easy Rise Walking Frame

Equipped with two-level hand grips, the Easy Rise Walker is designed to assist with the transition from sitting to standing.

- Made from sturdy yet lightweight aluminum tubing for easy maneuverability
- Non-slip, non-marking grips on front and rear legs
- One-touch folding mechanism ideal for storage and travel

ITEM# WA003

Gecko Starr / Reziprok Walker

With height-adjustable legs, supportive aluminum frame and an easy-to-fold, one-touch mechanism, the Gecko Starr Walker is the ideal space-saving mobility aid.

- Static, stable aluminium walker
- Ideal for safe walking
- Rubber tips prevent slipping
- Possible to use over the toilet basin

ITEM# 700100100 (Starr)

ITEM# 700100200 (Reziprok)



ITEM# WA009



ITEM# WA008

Item #	OVERALL DIMENSIONS CM (IN)				WEIGHT KG (LBS)		UOM	Color
	Overall Width Open	Overall Width Closed	Overall Length	Overall Height	Product Weight	Weight Capacity		
CLEVER-LITE WALKER - 5" CASTERS								
10230	65 (26)		65 (26)	80-100 (31-39)	9,1 (20.1)	136 (300)	1/cs	Black
FOLDING WALKING FRAME								
WA008	56 (22)		43 (17)	81-99 (32-39)	3 (6.6)	180 (397)	1/cs	Silver
WA009	64 (25)		46 (18)	88-98 (31-38,5)	3 (6.6)	180 (397)	1/cs	Silver

Item #	OVERALL DIMENSIONS CM (IN)				WEIGHT KG (LBS)		UOM	Color
	Overall Width Open	Overall Width Closed	Overall Length	Overall Height	Product Weight	Weight Capacity		
EASY RISE WALKING FRAME								
WA003	62 (24)		70 (28)		7,2 (15.9)	136 (300)	1/cs	Silver
GECKO STARR / REZIPROK WALKER								
700100100	61 (24)	12 (5)	47 (19)	81-98,5 (32-39)	2,5 (5.5)	120 (265)	1/cs	Silver
700100200	61 (24)	12 (5)	51 (20)	80-97,5 (31-38)	2,8 (6.2)	100 (220)	1/cs	Silver

Versatile and comfortable



Easy to unfold



ITEM# 10365-ADJ



ITEM# 10370BC

Cane with integrated seat

Capable of opening into a comfortable seat in one single motion, the Cane with an integrated seat is ideal for users with limited mobility.

- The legs create stability even on uneven surfaces and are height adjustable and stable.
- Height adjustable between 86 – 97 cm (34" - 38")
- Folds easily - 87 x 26 x 5 cm (34" x 10" x 2")

ITEM# 10365-ADJ



Gel Grip Cane

The soft, anatomical shaped fit of the Gel Grip cane provides optimal comfort while walking indoors and outdoors. This cane folds compactly for easy storage and travel.

- Optimal adjustment to body size
- Anti-slip rubber tip for added security
- Comfortable wrist strap
- Technical data and part number may differ depending on country

ITEM #	PRODUCT	COLOR	DESIGN
10365-ADJ	Cane w integr. seat	Champagne	
10370RC	Gel Grip Cane	Red Crackle	
10370BC	Gel Grip Cane	Blue Crackle	
10370BK	Gel Grip Cane	Black	

Specifications are available on page 33



ITEM# HCANE-BK-C2

Pivoting head



Folded



Handle



Feet



HurryCane® Freedom Edition™

Featuring eight convenient height adjustments, a 360 degree pivoting head and a free-standing design, the HurryCane Freedom Edition can accommodate any user.

- Anodized aluminum shaft
- Wrist strap and travel bag included
- Folds easily for convenient storage and travel
- Floor Display contains: 6 Black, 2 Red, 2 Blue, 2 Purple, 4 Heavenly Handles, 4 Packs of replacement feet
- Counter Display contains: 3 Red, 3 Purple, 3 Black, 3 Blue

ITEM# see table below



ITEM# HC-POPVR-E2-C2



ITEM# HC-POPVR-CA-C1

ITEM #	PRODUCT	COLOR
HCANE-BK-C2	HurryCane	Black
HCANE-BL-C2	HurryCane	Blue
HCANE-PR-C2	HurryCane	Purple
HCANE-RD-C2	HurryCane	Red
HC-HAND	HurryCane, Heavenly Handle	
HC-FEET	HurryCane, Replacement Feet	
HC-POPVR-CA-C1	Floor display	
HC-POPVR-E2-C2	Counter display	

Specifications are available on page 33



ITEM# WS016

ITEM# WS016L



ITEM# 10400-30

ITEM# 10452-30

ITEM# 10483

Adjustable legs



ITEM #	COLOR	SIZE	DESIGN
WS016	Silver	M	
WS016L	Silver	L	
10400-30	Silver	Adult	
10452-30	Silver		
10481	Red		
10482	Blue		
10483	Silver		

Specifications are available on page 33

Quad Cane M / L

The lightweight and comfortable Medium Quad Cane has a four-prong foot for added stability, and an offset ergonomically designed handle that evenly distributes the user's weight while walking.

- Height adjustable
- Secure and comfortable grip
- Strong aluminium shaft

ITEM# WS016 / WS016L

Aluminum Crutches

Featuring a double extruded center tube for additional weight-bearing strength, the Lightweight Aluminum Crutches provide optimal durability and stability.

- Non-skid, jumbo-size vinyl contoured tips provide excellent traction
- Easy wing nut adjustments for hand grips
- Latex-free, comfortable, durable accessories
- Easy pushpin height adjustments in 2.5 cm (1") increments
- Weight capacity of 159 kg (350 lbs.) on aluminum crutch-differs for pediatric

ITEM# 10400-30

Forearm Crutches Double Adjustable

The Double Adjustable Forearm Crutches are slip resistant to provide the user with optimal safety and stability, and the double height adjustment of the legs and cuffs are ideal for recovery.

- Crutch height adjusts in increments of 2.5 cm (1") for optimal sizing
- Slip resistant, gray rubber grips provide stability
- Lightweight and latex-free

ITEM# 10452-30

Forearm Crutches Adjustable

Equipped with adjustable legs, these forearm crutches accommodate users of all shapes and sizes.

- Lightweight aluminium
- Latex free
- Anti-slip grip for added security

ITEM# 10481 / 10482 / 10483



Item #	OVERALL DIMENSIONS CM (IN)		SEAT CM (IN)			WEIGHT KG (LBS)			UOM	Color
	Overall Height		Seat Width	Seat Depth	Seat Height	Product Weight	Carton Shipping Weight	Weight Capacity		
GEL GRIP CANE										
10370RC	81-93 (32-37)					0,4 (.8)		130 (287)	1/cs	Red Crackle
10370BC	81-93 (32-37)					0,4 (.8)		130 (287)	1/cs	Blue Crackle
10370BK	81-93 (32-37)					0,4 (.8)		130 (287)	1/cs	Black
CANE WITH INTEGRATED SEAT										
10365-ADJ			21 (8)	21 (8)	48-56 (19-22)	0,4 (.8)		100 (220)	1/cs	Champagne
HURRYCANE® FREEDOM EDITION™										
HCANE-BK-C2	77-95 (30.5-37.5)							175 (350)	6 pcs	Black
HCANE-BL-C2	77-95 (30.5-37.5)							175 (350)	6 pcs	Blue
HCANE-PR-C2	77-95 (30.5-37.5)							175 (350)	6 pcs	Purple
HCANE-RD-C2	77-95 (30.5-37.5)							175 (350)	6 pcs	Red
HC-HAND									1/cs	
HC-FEET									3cs	
HC-POPVR-CA-C1									1/cs	Silver
HC-POPVR-E2-C2									1/cs	Silver
QUAD CANE M										
WS016	66-92 (26-36)					0,75 (1.6)		115 (253.5)	1/cs	Silver
QUAD CANE L										
WS016L	79-104 (31-41)					1,3 (2.8)		115 (253.5)	1/cs	Silver
ALUMINUM CRUTCHES										
10400-30	114-135 (45-53)							159 (351)	Pair	Silver
FOREARM CRUTCHES DOUBLE ADJUSTABLE										
10452-30	65-91 (26-36)							127 (280)	Pair	Silver
FOREARM CRUTCHES ADJUSTABLE										
10481	84-99 (33-39)					1,6 (3.5)		100 (220)	1/cs	Red
10482	84-99 (33-39)					1,6 (3.5)		100 (220)	1/cs	Blue
10483	84-99 (33-39)					1,6 (3.5)		100 (220)	1/cs	Silver



▶ **Patient Room**

Beds

Mattresses

Pressure Area Care

Plinths

Patient Lifts

Riser Recliner

Overbed Tables

Daily Aids



Casa Med Classic FS with Side Rails

Offering a classic, attractive and non-clinical finish, the Casa Med Classic FS with Side Rails features a Dewert motor, 8-button handset and controls. The electric backrest, leg rest, knee brake and lower leg elevation make it easy to customize to the needs of the caregiver and user. Available in a variety of colors, this standard height bed (40-80 cm/16"-32") enables caregivers to care for users successfully without bending and causing back injuries.

- Emergency lowering by quick action bolt mechanism
- Four brake casters (100mm/3.9")
- Four section mattress support (mesh or wooden slat)
- Quick-release actuator points
- Extra-high side rails available
- Tool-free assembly and dismantle

ITEM #	POWER CORD	MATTRESS SUPPORT
19849SED3BE	UK	Metal Mesh
19847SED3BE	UK	Wooden Slat
19849SEUROED3BE	EU	Metal Mesh
21847ED3BE	EU	Wooden Slat

Specifications are available on page 47



Casa Med II with Side Rails

The Casa Med II Bed is elegantly finished with fully enclosed motors and real wood head and foot boards. It offers an additional scissor-lift raising mechanism to provide greater lifting stability and an increased safe working load. With a Dewert motor, 8-button handset and controls, this bed provides an electric backrest, leg rest, knee brake and lower leg elevation.

- Emergency lowering by quick action bolt mechanism
- Four brake casters (100mm/3.9")
- Four section mattress support (mesh or wooden slat)
- Quick-release actuator points
- Extra-high side rails available
- Tool-free assembly and dismantle
- Transport system and lock

ITEM #	POWER CORD	MATTRESS SUPPORT
21871ED3BE	UK	Metal Mesh
19865ED3BE	UK	Wooden Slat
21876ED3BE	EU	Wooden Slat

Specifications are available on page 47



Casa Med Ultra FS with Side Rails

The Casa Med Ultra FS is a competitively priced entry-level bed with all the high levels of functionality you would expect from a Casa Med bed. With a Dewert motor, 8-button handset and controls, this bed provides an electric backrest, leg rest, knee brake and lower leg elevation.

- Emergency lowering by quick action bolt mechanism
- Four brake casters (100mm/3.9")
- Four section mattress support (mesh or wooden slat)
- Quick-release actuator points
- Extra-high side rails available
- Tool-free assembly and dismantle
- Transport system and lock

ITEM #	POWER CORD	MATTRESS SUPPORT
19909SED3BE	UK	Metal Mesh
19908SED3BE	UK	Wooden Slat
19909SEUROED3BE	EU	Metal Mesh
19908SEUROED3BE	EU	Wooden Slat

Specifications are available on page 47



Elitis+ Standard Medical Bed

Available in three models, the Elitis+ Bed is a clinical bed assuring simplicity, features and stability. Its mattress support is height adjustable from 36-83 cm (14"-33"). The optional numerous accessories can help to satisfy most patient's needs, and a choice of finishes helps to personalize the environment.

Three models:

Two electrical functions:

- Electrical height adjustment
- Electrical backrest
- Manual reclining legrest with racks (2)

Three electrical functions:

- Electrical height adjustment
- Electrical backrest
- Electrical reclining legrest

Three electrical functions with kneebreak:

- Electrical height adjustment
- Electrical backrest
- Manual reclining legrest with electrical kneebreak and adjustment of the foot section with racks

ITEM #	
S6100400	

Specifications are available on page 47



INNOV8 iQ Hospital Bed

The Innov8 iQ high performance hospital bed has successfully brought together an in-depth understanding of the changing requirements of the acute sector with state of the art manufacturing and assembly methods. This innovative ward bed achieves the optimum combination of safety, performance, functionality and cost efficiency.

- Low height reduces risks associated with patient falls and increases patient independence.
- Backrest and mattress angle indicators
- Removable head and foot ends available in a number of color options
- Movable catheter hooks
- Integral bed extension combines ease of use and functionality



ITEM #	SIDE RAILS
IQ/AQ/SR	With Standard Side Rails
IQ/AQ/SSR	With Split Side Rails

Specifications are available on page 47



ITEM #	SIDE RAILS
INNOV8/LOW/AQ/LINK/SR	With Standard Side Rails

Specifications are available on page 47

INNOV8 Low Bed

Designed for safe access and exit, the Innov8 Low Bed features a minimum platform height of just 218 mm (8.6") – one of the lowest available on the market today. Designed to provide total flexibility in terms of bed specification for applications from hospital beds through to high-dependency environments.

- Exceptional low-to-high range- 22 to 70 cms (8.7" to 27.6")
- Smooth, electrically operated backrest, height and kneebreak (with battery backup)
- Trendelenburg, reverse Trendelenburg, auto-contour and auto-regression plus auto-cardiac chair position
- Simple-to-use patient handset and nurse control panel with handset lockout
- Dual-sided CPR levers
- Extendable mattress platform length for taller patients





Simple Air Mattress

The Simple Air Mattress features an easy-to-use control panel with a range of comfort settings to accommodate varying comfort needs. The easy-access and simple-to-use CPR tag allows rapid deflation of the mattress within 10 seconds to aid emergency procedures.

- Adjustable weight setting
- Adjustable dynamic cycle times (10,15,20,25 min)
- Adhesive pads
- Audible and visual alarm
- Autofirm setting
- External conventional fuse
- Full-length hose sleeve
- 360° zip



ITEM #	
1001270001	
Specifications are available on page 47	



Memaflex Advance Mattress

The Memaflex Advance mattress features a high-density visco-elastic foam surface, ideal for delivering exceptional patient comfort and a high degree of pressure reduction. The temperature-sensitive foam conforms to the shape of the body, creating a greater contact area and providing effective pressure redistribution.

- Internal air channels in the foam's surface help to reduce the build up of heat and moisture on the skin
- Smooth multi-stretch cover allows the patient to glide over the surface of the mattress
- Individually castellated cells conform independently, allowing the mattress to profile effectively

ITEM #	
2001020001	
Specifications are available on page 47	



Hybrid-Power 2-in-1 Mattress

The Hybrid-Power 2-in-1 Mattress combines the pressure redistribution properties of high-specification foam with alternating pressure relief. The visco-elastic foam conforms to the shape of the body, while the alternating cells assist with pressure relief.

- Cover can withstand a higher concentration of cleaning substances and offers greater resistance to abrasion
- High-frequency welded seams protect against fluid ingress
- 360° zip to prolong longevity
- Automatic lockout feature
- Low pressure alarm and indicator
- Compact and lightweight power unit
- Power failure alarm and indicator
- Non-slip base – prevents the mattress from sliding out of position, giving reassurance to the patient and caregiver



ITEM #	
1001290001	
Specifications are available on page 47	

Permafex Plus Advance Mattress

The Permafex Plus Advance mattress features a castellated support surface within a U-shaped foam base. This unique design provides maximum support during patient transfers and delivers exceptional pressure reduction.

- Waterproof, vapour-permeable cover can be easily wiped clean
- High-frequency welded seams reduce the risk of fluid ingress
- Zip is protected by a double sided, full length welded flap
- Cover can be removed and laundered at 80°C




ITEM #	
2001030003	
Specifications are available on page 47	



Community Advance Mattress

Ideal for use in a nursing home, the Community Advance is an ideal entry-level mattress for patients at some risk of developing pressure ulcers, who require a high level of comfort.

- Cover and foam can be replaced individually
- Waterproof, vapour-permeable cover can be easily wiped clean
- Concealed zip helps to reduce the risk of fluid ingress
- Cover can be removed and laundered at 80°C
- Also available in a double size mattress



ITEM #	
2001090001	
Specifications are available on page 47	

Acclaim Flow Non Powered Hybrid Mattress

The Acclaim Flow is a non-powered mattress replacement system using Sidhil's Air Flow Technology, a clever combination of air cells and foam developed to achieve optimum pressure reduction. The system features automatically adjusting valves across 10 individual latitudinal air cells which respond to the patient's weight and movements without the need for a powered pump.

- Zoned, castellated foam overlay – individual heel, head and body zones facilitate appropriate localized pressure reduction
- Firm edges for ease of patient ingress/egress and transfers
- Heavy-duty, anti-slip base, vapour permeable and waterproof multi-stretch cover with welded seams and waterfall flap
- Internal cells incorporated securely into foam



ITEM #	
MAT/ACCL/FLOW	
Specifications are available on page 47	



Artemis Dynamic Mattress

This dual therapy system offers both alternating therapy for effective pressure relief and constant low pressure to reduce peak pressures. The combination of therapies gives the Artemis optimum flexibility to meet diverse patient requirements within ITU, ward or community care environments.

- Two-cell alternating dynamic pressure therapy, full mattress replacement system
- Fully automatic pressure adjustment – instantly reacts to changes determined by patient weight and positioning
- Zoned mattress system for optimum outcomes - static head cells and narrow heel cells
- Figure 8 cell construction with the sub-cell remaining inflated to reduce risk of the patient bottoming out
- Patient egress alarm – alerts caregiver when no patient weight is detected on the mattress
- Easy-to-use digital control system incorporates alternating, constant low pressure, pulsation and max inflate modes
- Dual compressor for quick inflation / Rapid deflate CPR

ITEM #	
DYN/DIG/ARTEMIS	
Specifications are available on page 47	



MAP Continuous Bedside Pressure Mapping System

M.A.P™ is the very first Continuous Bedside Pressure Mapping System, which can be used on any mattress system to provide 24/7 data on pressure levels developing between patient and support surface.

Presented graphically on a screen, this real time information on pressure distribution gives nursing staff and carers the accurate detail they need to reposition patients effectively, helping to reduce the incidence of pressure ulcers.



ITEM #	
MAP/100/S	
Specifications are available on page 47	



Serenade Gel Cushion with Straps

The Serenade Gel Cushion is a silicone gel flotation system developed to reduce pressure for high-risk users.

- Strap to affix to chair, secures firmly in place
- Maximum weight capacity 160kg (353 lbs.)

ITEM #	COLOR
MAT/SERE/GEL/STRAP	Midnight Blue
<i>Specifications are available on page 47</i>	



Memaflex Cushion

The Memaflex Cushion features a high-density, temperature-sensitive, visco-elastic cushion, allowing it to conform to the user's shape and maximizing body contact area. It was designed for patients at high risk of developing pressure ulcers and for those patients who require comfort as well as pressure prevention.

- Waterproof, multi-stretch, vapor-permeable polyurethane cover
- High-frequency welded seams and fully enclosed zip
- Non-slip base

ITEM #	COLOR
2002140001	Midnight Blue
<i>Specifications are available on page 47</i>	





Salisbury Variable Height Couch

The Salisbury Variable-Height Couch is an electrically controlled, multi-function couch with four-braked, swivel casters for ease of maneuverability and safety.

- Electric handset control
- Swing away side rails and paper roll holder fitted as standard
- Adjustable head and backrest section

ITEM #	COLOUR
COU32/(Colour)/1	see colour scheme
<i>Specifications are available on page 47</i>	

Available Colours



Doherty Two-Section Plinth

Developed specifically for applications within GP practices, podiatry and physiotherapy clinics, this electrically operated, two-section plinth is a high-specification unit. Designed for optimum flexibility, it offers both head-down tilt and feet-up positions. Delivered fully assembled, a breathing hole comes standard on the headrest section.

- Height range from 46 cm to 85 cm (18" to 33.5")
- Easi Lift, gas-assisted backrest
- Electric handset control for height adjustment
- Twin retractable casters for ease of movement
- Distinctive Doherty style upholstery with rounded corners for hygiene and cleanliness
- Dual-sided operation for positioning flexibility



ITEM #	COLOUR
PLE01/(Colour)/1	see colour scheme
<i>Specifications are available on page 47</i>	

Item #	OVERALL DIMENSIONS CM (IN)			PUMP				WEIGHT KG (LBS)				UOM
	Width	Length	Height	Power (V, Hz, A)	Airflow (LPM)	Cycle Time (Minutes)	Type	Product Weight	Carton Shipping	Weight Capacity	Heaviest compo-	
CASA MED CLASSIC FS WITH SIDE RAILS												
All item#	102 (40)	213 (84)	40-80 (16-32)	230 50/60Hz				91.5 (201.7)		140 (309)	20 (44)	1/cs
CASA MED ULTRA FS WITH SIDE RAILS												
All item#	213 (84)	213 (84)	40-80 (16-32)	230 50/60Hz				91.5 (201.7)		140 (309)	20 (44)	1/cs
CASA MED II WITH SIDE RAILS												
All item#	101 (40)	206 (81)	40-80 (16-32)	230 50/60Hz				149 (328)		165 (364)	40 (88)	1/cs
ELITIS+ BED												
S61000400	91 (36)	205 (81)	36-83 (14-33)					80 (176.4)		170 (375)		1/cs
INNOV8 IQ HOSPITAL BED												
IQ/AQ/SR & IQ/AQ/SSR	103 (41)	228 (90)	34-75 (13-30)	230 50/60Hz				130/146		255 (562)		1/cs
INNOV8 LOW BED												
INNOV8/LOW/AQ/LINK/SR	102 (40)	218 (86)	22-71 (9-28)	230 50/60Hz				130 (286)		255 (562)		1/cs
SIMPLE AIR MATTRESS												
1001270001	90 (35)	200 (78.7)	20 (7.8)	220, 50/60	8	10/15/20/25	Compressor	7 (15.4)		180 (396)		1/cs
HYBRID-POWER 2-IN-1 MATTRESS												
1001290001	88 (34.6)	198 (78)	16 (6.3)	230, 50, 0.5	8	10	Compressor	16 (35)		254 (560)		1/cs
MEMAFLEX ADVANCE MATTRESS												
2001020001	88 (34.6)	198 (78)	16 (6.3)					16 (35.3)		247.7 (546)		1/cs
PERMAFLEX PLUS ADVANCE MATTRESS												
2001030003	88 (34.6)	198 (78)	16 (6.3)					14 (30.8)		247.7 (546)		1/cs
COMMUNITY ADVANCE MATTRESS												
2001090001	88 (34.6)	198 (78)	12.5 (5)					12 (26.5)		165 (363)		1/cs
ARTEMIS DYNAMIC MATTRESS												
DYN/DIG/ARTEMIS	88 (34.6)	198 (78)	20 (7.8)	230 50Hz		10	Compressor	15 (33)		248 (546)	11	1/cs
ACCLAIM FLOW NON POWERED HYBRID MATTRESS												
MAT/ACCL/FLOW	88 (34.6)	199 (78)	15.2 (6)				Self Regulating	11 (24.3)		230 (507)		1/cs
MAP CONTINUOUS BEDSIDE PRESSURE MAPPING SYSTEM												
MAP/100/S												1/cs
SERENADE GEL CUSHION WITH STRAPS												
MAT/SERE/GEL/STRAP	43 (17)	43 (17)	7.5 (3)					2 (4.4)		160 (352)		1/cs
MEMAFEX CUSHION												
2002140001	43 (17)	43 (17)	5 (2)					3 (6.6)		127 (280)		1/cs
DOHERTY TWO-SECTION PLINTH												
PLE01/(Colour)/1	60 (24)	193 (76)	46-85 (17-33)	230 50Hz				78 (172)		228 (502)		1/cs
SALISBURY VARIABLE HEIGHT COUCH												
COU32/(Colour)/1	61 (24)	185 (73)	52-85 (20-33)	230 50Hz				70 (154)		160 (352)		1/cs



Hermes Lifter

With a 250 kg (550 lb.) weight capacity, the Hermes Lifter features a cross bar equipped with four secure points for reinforced stability. This lifter is also equipped with double-roller, 100 mm (4") wheels with solid tire casing for exceptional mobility.

- High-capacity, electrical moving base
- Stored vertically
- Low clearance for easy access underneath beds

ITEM #
SR2310000
<i>Specifications are available on page 53</i>

Samssoft 175 V2 Lifter

The electrically operated Samssoft 175 Lifter has manual pedals, enabling the chassis legs to be opened with ease and wheeled around a chair. The lightweight frame is easy to move into position, allowing the user to achieve a lower lifting position from the floor when transferring patients.

- Stored vertically
- No disassembly necessary for transport
- Low clearance for easy access underneath beds
- 175 kg (385 lbs.) weight capacity

ITEM #	SPECIFICATION
SR2810100EU	EU plug
SR2810100UK	UK plug
<i>Specifications are available on page 53</i>	



Casa 180i Lifter

The Casalift Hoist is a lightweight and versatile lifting device for use when transferring individuals.

- Easy manoeuvrability and access
- Lightweight frame is easy to move into position
- High lifting range and safe working load
- Ergonomic, user friendly handset



ITEM #
18001BE
<i>Specifications are available on page 53</i>



Chicago Riser Recliner

The Chicago Riser Recliner features a single, quiet motor to activate a distinct lift, rest and recline action with a simple-to-operate, two-button handset. With a stylish design and four color options to fit virtually any decor, this riser recliner is ideal for individuals who have difficulty standing up and sitting down due to restricted mobility.

- Ergonomically contoured seat and back
- Padded upholstery for added comfort and support
- Excellent snooze position with backrest reclining to 64° for comfort and relaxation

ITEM #	COLOR
CLR11BE	Beige
CLR11GO	Gold
CLR11GR	Green
CLR11T	Terracotta

Specifications are available on page 53



ITEM# CLR11GO



ITEM# CLR12BE

Nevada Riser Recliner

The Nevada Riser Recliner features two discreet motors that enable the legrest to be elevated and lowered as well as the backrest to be reclined independently to provide more flexible support options. Ideal for individuals who have difficulty standing up and sitting down due to restricted mobility, this riser recliner is also available in petite size for individuals 1.7m (5'6") tall or less.

- Includes simple-to-operate, four-button handset
- Available in luxurious Rembrandt fabric with four color options

ITEM #	COLOR
CLR12BE	Beige
CLR12GO	Gold
CLR12GR	Green
CLR12T	Terracotta

Specifications are available on page 53

Over Bed/Chair Table

With stable steel construction, the Mobile Open Toe Table offers height and tilt adjustment enabling the top to be set at varying heights and angles as needed.

- Durable plastic coated 2.5 cm (1") square tubing
- Easy-to-clean laminated teak top with raised edges

ITEM #	
SA4211000	

Specifications are available on page 53



ITEM# SA4211000



ITEM# 948

Deluxe Table Twin Top

The handy Deluxe Table Twin Top is ideal for use in the bedroom or living room. Featuring a small side table for drinks while the main table is tilted for reading, this overbed table has four casters that allow for easy maneuvering and positioning.

- Simple to assemble
- Melamine-coated top with a raised side edges to prevent items from falling
- Height and tilt adjustment enabling top to be set at varying

ITEM #	
948	

Specifications are available on page 53



Ring Pull Can Opener

- Assists individuals with reduced grip to open cans with less effort
- Easy-to-grip rubber handle
- Suitable for left or right handed users

ITEM #	MEKD002BLK-WHT
--------	----------------



Universal Cup

- Suitable for individuals with reduced grip or limited sensation in their hands
- Two large side handles allow the user to use whole hand to grip handle or to put hand through handle to hold the cup, enabling more stable grip
- Spout and anti-splash lid prevent spillage

ITEM #	RTL3515
--------	---------



Pedal Trainer

- Enables individuals to exercise legs and arms while sitting comfortably
- Easy-to-use while watching TV or reading

ITEM #	10273KDR
--------	----------



Hip Kit

- Long Handled Bath Sponge
- Poly Stocking Aid
- Hand-held Reacher 81 cm (32")
- Plastic Shoe Horn
- Black Elastic Shoe Laces

ITEM #	RTL9506-7
--------	-----------



Reacher

- Enables individuals to reach high shelves safely
- Special rubber lips enable grabbing even the smallest objects
- Robust construction yet light in the hand

ITEM #	PC23583
--------	---------



Handled Sponges

- 61 cm (24") strong, long, acrylic handle
- 76 cm (30") strong, extra long, acrylic handle

ITEM #	HANDLED SPONGES
RTL1020	Long
RTL1027	Extra Long



Swivel Cushion

- Cushion rotates 360° to support movement
- Suitable for all common seats
- Diameter: 40 cm (15.7")
- 120 kg (265 lb.) weight capacity

ITEM #	SWIV001R
--------	----------



Eating Utensils

- Utensils suited for individuals with arthritis, weak or limited handgrip
- Ergonomically designed contoured grip allows for comfortable fit in almost any hand while adding to stability

ITEM #	UTENSIL
RTL1410	Fork
RTL1411	Spoon
RTL1413	Knife

Item #	HERMES LIFTER			SAMSOF 175 V2 LIFTER			CASO 180i Lifter			CHICAGO RISER RECLINER			NEVADA RISER RECLINER			OVER BED/CHAIR TABLE			DELUXE TABLE TWIN TOP		
	Width Open	Width Closed	Length	Height	Width	Depth	Height	Product Weight	Carton Shipping Weight	Weight Capacity	UOM	Ships from	Width	Depth	Height	Product Weight	Carton Shipping Weight	Weight Capacity	UOM	Ships from	
SR2310000	74 (29)		136 (54)	200 (79)				55,0 (121.3)		250 (551)	1/cs										
SR2810100EU	61 (24)		128 (50)	159 (63)				36,0 (79.4)		175 (386)	1/cs	FR									
SR2810100UK	61 (24)		128 (50)	159 (63)				36,0 (79.4)		175 (386)	1/cs	FR									
SR2910000	61 (24)		94 (37)	88 (35)				39,0 (86.0)		175 (386)											
CLR11BE	83 (33)		93 (37)	108 (43)	50 (20)	52 (20)	48 (19)			150 (331)											
CLR11GO	83 (33)		93 (37)	108 (43)	50 (20)	52 (20)	48 (19)			150 (331)											
CLR11GR	83 (33)		93 (37)	108 (43)	50 (20)	52 (20)	48 (19)			150 (331)											
CLR11T	83 (33)		93 (37)	108 (43)	50 (20)	52 (20)	48 (19)			150 (331)											
CLR12BE	86 (34)		93 (37)	110 (44)	51 (20)	54 (21)	48 (19)			150 (331)											
CLR12GO	86 (34)		93 (37)	110 (44)	51 (20)	54 (21)	48 (19)			150 (331)											
CLR12GR	86 (34)		93 (37)	110 (44)	51 (20)	54 (21)	48 (19)			150 (331)											
CLR12T	86 (34)		93 (37)	110 (44)	51 (20)	54 (21)	48 (19)			150 (331)											
SA4211000					40 (16)	45 (18)															



▶ **Bath Safety**

Wheeled Commodes

Commodes

Raised Toilet Seats

Bath Lifts

Bath Benches & Steps

Bath Seats & Stools

Shower Chairs

Shower Wheelchairs

TRS 130 Wheeled Commode

Comfortable, safe and durable, the TRS 130 Wheeled Commode features a removable seat and cushion for access to the receptacle. This commode can be used with the receptacle or easily positioned above a toilet basin for use without.

- Fold-down armrest simplifies transfer
- Stable, retractable and removable footrests
- Seat height adjustment through telescopic legs
- Rear caster parking brakes for safety
- Simple, tool-free assembly
- Easy to clean

ITEM #	SEAT WIDTH CM (IN)
530100000	44 (17)
<i>Specifications are available on page 68</i>	



Deluxe Drop Arm Commode

With a 300 kg (661 lb) weight capacity, the Deluxe Drop Arm Commode features an easy-to-release arm mechanism, which allows for safe lateral user transfers to and from commode. In addition, the large, durable snap-on seat can be removed for easy cleaning, and the covered, 10-liter bucket helps to prevent spills and confine odors.

- Large backrest for added user comfort and support
- Durable, heavy duty, white powder coated steel tubing
- Padded armrests for additional comfort and support
- Bucket easily slides in and out for emptying and cleaning

ITEM #	SEAT WIDTH CM (IN)
11135-1	58 (23)
<i>Specifications are available on page 68</i>	



TS 130 Commode

The comfortable and durable TS 130 Commode features a large, soft seat and backrest to provide added support. The telescopic legs enable seat height adjustment, and the rubber feet keep the commode stable while in use.

- Seat and cover can be lifted and removed
- Standard receptacle
- Easy to clean

ITEM #	SEAT WIDTH CM (IN)
530100100	43 (17)
<i>Specifications are available on page 68</i>	




Bariatric Folding Commode

With a 295 kg (650 lb) weight capacity, the Bariatric Folding Commode easily folds for storage and transportation, and assembles tool free in three easy steps.

- Durable, plastic 12-quart bucket with snap-on seat and lid
- Plastic armrests for added comfort
- Removable tool-free back

ITEM #	SEAT WIDTH CM (IN)
11117N-2	35 (14)
<i>Specifications are available on page 68</i>	



Ticco 2G Raised Toilet Seat

Available with or without a lid, the lightweight Ticco 2G Raised Toilet Seat features a large opening with cutaways for enhanced hygiene.

- Simple, tool-free assembly
- Fixing screws for precise hold
- Smooth surfaces for simple cleaning



ITEM #	SEAT WIDTH CM (IN)	COVER
523400000	39 (15)	✓
523300000	39 (15)	
<i>Specifications are available on page 68</i>		




Elevated Toilet Seat

Designed to fit on existing toilet bowls, the Elevated Toilet Seat features removable arms to facilitate side transfers and aids in sitting down and standing up from the toilet.

- Heavy-duty molded plastic construction for additional strength and durability
- Foam padding across the entire length of the armrests for added comfort
- Improved locking device to ensure a safe, secure fit
- Front seat recess for personal hygiene access
- Tool-free attachment and non-skid, anti-slip for user safety




ITEM #	SEAT WIDTH CM (IN)
12027RA	43 (17)
<i>Specifications are available on page 68</i>	



TSE 150 Raised Toilet Seat

The TSE 150 Raised Toilet Seat features two hygienic cutaways and an integrated seat. Retractable, slip-resistant handle grips provide stability and security when mounting and dismounting the seat.

- Recessed grip on the cover makes it easy to open
- Simple, tool-free assembly
- Three-position height adjustment
- Smooth surfaces for simple cleaning



ITEM #	WIDTH CM (IN)	ARMS
5201000	56 (22)	✓
5200500	46 (18)	
<i>Specifications are available on page 68</i>		



Cosby Toilet Frame

The free-standing, adjustable Cosby Toilet Frame surrounds the toilet and provides support when getting on and off toilet.

- Available with molded plastic armrests or padded polyurethane armrests for a soft, comfortable grip
- Sturdy steel frame coated with Duradip® tough, peel-resistant plastic to combat corrosion and improve longevity
- Non-slip, non-marking feet ensure stability with a floor fixed version also available for additional safety
- Stackable to maximise storage space
- Suitable for thermal decontamination in automated washing machines

ITEM #	WIDTH CM (IN)
157	48-61 (19-24)
<i>Specifications are available on page 68</i>	



Bellavita 2G Bath Lift

Weighing 21 pounds, the Bellavita 2G is one of the lightest bath lifts on the market. Featuring an ergonomic seat and backrest, this lift provides comfort and safety to the user while bathing.

- Wide side flaps available as accessories
- Quick-release front suction cups
- Tested for electromagnetic compatibility (EMC)
- No height adapter required for baths up to a depth of 48 cm (19")
- Five-year warranty on bath lift, two-year warranty on



ITEM #	PRODUCT
475100312	With classic blue cover
475150312	With comfort blue cover

Specifications are available on page 68

Bellavita 2G Accessories



Headrest

This padded headrest is height adjustable for additional comfort.
ITEM# 410200312 Blue



Swivel and Transfer Aide Vitaturn L

This swivel transfer aide facilitates with transfers and positioning to accommodate the user's individual needs.
ITEM# 412100200



Side Flap Deflectors

These deflectors mount easily and prevent the side flaps from getting caught in the handles.
ITEM# 400570



Belt System Two-Point Lap Belt

The belt system offers safety and comfort while bathing and installs easily.
ITEM# 460900100



Foam Cleaner

This foaming cleaner is specifically for the plastic surfaces of the Bellavita 2G.
ITEM# 410900100



Bellavita Classic Blue Cover Set

This non-slip cover set provides comfortable seating while bathing.
ITEM# 460800312



Suction Cup Set

The suction cup set holds fast when needed and tab vents for release upon bending.
ITEM# 460900225 Clear



Bellavita Comfort Blue Cover Set

Featuring a soft, non-slip material, this set includes covers for the side flaps that are washable and easy to clean.
ITEM# 460900312



Widening Side Flaps

These side flaps widen the seating and transfer surface by 5, 7.5 or 10 cm (2", 3" or 4").
ITEM# 460900227 SIZE 5 cm
ITEM# 460900228 SIZE 7.5 cm
ITEM# 460900229 SIZE 10 cm



Seat Shortener

This Seat Shortener shortens the seating surface for secure positioning.
ITEM# 460950252



Swivel Cushion Vitaturn K

With a height of 12 mm (.5"), the machine washable swivel seat cushion accessory assists individuals getting in and out of the bathtub.
ITEM# 412100100



Belt System Four-Point Chest Belt

Featuring an easy-to-install design, the harness offers additional safety and comfort to the user while bathing.
ITEM# 460900150



Benny Bath Bench

Featuring a grab bar and slip-resistant seat surface, the Benny bath bench is ideal for safe and easy bathtub entry.

- Equipped with drainage holes in the seat
- Adjusts easily to fit any tub width
- Includes rubber feet for added safety and protection from scratching the tub
- Tool-free installation

No-scratch rubber foot



ITEM #	LENGTH CM (IN)
550100100	69 (27)
550100200	74.5 (29)

Specifications are available on page 68

Sky Bath Step

Equipped with a large, slip-resistant step surface, the Sky bath step simplifies and provides additional safety when getting in and out of the bathtub.

- Stacks easily to raise the step height
- Easy-to-clean
- Weighs less than a pound

Easily stackable



Drainage holes



Individual adjustment



ITEM #	WIDTH CM (IN)
RTL12068	50 (20)

Specifications are available on page 68



Suction Cup 12" Grab Bar with Indicator

Equipped with a color indicator that shows when there is a safe suction and when the suction needs to be readjusted, the 12" grab bar is ideal for bath safety. The release levers make installation and removal quick and easy.

ITEM# 13082-3



BEN Swiveling Bath Seat

The fully-assembled BEN swivel bath seat can be positioned in 90-degree intervals for added comfort and safety while bathing.

- Equipped with drainage holes and an automatic latching system for added security
- Slip-resistant rubber mounting arms
- Capable of adjusting up to 59 cm (23") to accommodate varying bathtub widths

Adjustable screw with cap



Slip-resistant rubber ends



90 degree rotation



ITEM #	SEAT WIDTH CM (IN)
550600100	44 (17)

Specifications are available on page 68

Dino Shower Stool

Capable of adjusting to eight seat heights, the compact Dino shower stool offers a safe and stable bathing aide option

- Lightweight and does not require assembly
- Structured surface for secure hold
- Stable thanks to rubber feet

Adjustable legs



ITEM #	SEAT DIAMETER CM (IN)
550500100	32 (13)

Specifications are available on page 68



DKS 130 Folding Shower Chair

Featuring a sturdy wall-mounted and folding design, the DKS 130 shower chair provides safety and stability while showering.

- Adjustable feet for optimal stability
- Large seat with drainage holes
- Seat adjusts to four different heights
- Equipped with rust-resistant aluminum legs



ITEM #	SEAT WIDTH CM (IN)
550700100	41 (16)
<i>Specifications are available on page 68</i>	

Deluxe Shower Chair

The Deluxe shower chair features a wide seat to accommodate individuals of all shapes and sizes and provide additional safety and stability while bathing.

- Equipped with a leg crossbar for added stability
- Large seat with drainage holes
- Seat can adjust to five different heights



ITEM #	SEAT WIDTH CM (IN)
12009-2	68.6 (27)
<i>Specifications are available on page 68</i>	

DS 130 Shower Chair

The DS shower chair provides additional comfort, safety and stability while showering.

- Seat can adjust to five heights
- Equipped with a hygienic cutaway
- Large seat with drainage holes
- Armrest available separately



ITEM #	SEAT WIDTH CM (IN)	ARMREST
550200200	41 (16)	✓
<i>Specifications are available on page 68</i>		



Foxton Stationary Shower Chair

This free-standing shower chair is equipped with a backrest and armrest for additional support while bathing.

- Polypropylene molded seat with drainage holes
- Equipped with an H section base for added strength
- Easy e-clip height adjustment on some models for optimal sizing
- Plastic coated with rust-resistant Duradip

ITEM #	SEAT WIDTH CM (IN)
140	44 (17)
<i>Specifications are available on page 68</i>	

Uppingham Shower Commode

The Uppingham shower commode allows the caregiver to wheel the individual to and from the toilet or shower without the need for transfers.

- Padded, moisture-repellent polypropylene molded seat, backrest and armrests
- Plastic coated with peel-resistant Duradip
- Equipped with a commode bucket



ITEM #	SEAT WIDTH CM (IN)
172	42 (17)

Specifications are available on page 68

DuoMotion Shower Wheelchair

The DuoMotion shower wheelchair offers adjustable legrests and armrests and features 61 cm (24") wheels.

- Features a hygienic front cutaway and cover
- Half axles and brakes available separately
- Seat adjusts to four positions
- Equipped with standard receptacle
- Machine washable, soft backrest



ITEM #	PRODUCT	SEAT WIDTH CM (IN)	VERSION
530500100	DuoMotion 24	44 (17)	Self Propel
530500200	DuoMotion	44 (17)	Transit

Specifications are available on page 68

24" drive wheels Adjustable footrests



Aston Shower Commode

The Aston shower commode provides assistance when transferring an individual to and from the toilet or shower. This shower commode also features a self-propelling rear kit that can be attached to enable independent use.

- Can be used solely as a commode
- Padded, moisture-repellent polypropylene molded seat, backrest and armrests
- Features a hygienic front cutaway
- Lightweight, rust-proof aluminum frame for easy maneuvering and handling



Transit Version



ITEM #	PRODUCT	SEAT WIDTH CM (IN)
CWC003	Aston	46 (18)
CWC003HD	Aston Heavy Duty	56 (22)

Specifications are available on page 68



DuoMotion XL Shower Wheelchair

The DuoMotion XL shower wheelchair is equipped with an additional crossbar for added stability and security.

- Features a hygienic front cutaway and cover
- Seat adjusts to four positions
- Equipped with standard receptacle
- Machine washable, soft backrest
- Available with folding armrests and adjustable legrests that are equipped with heel straps
- Includes 61 cm (24") wheels with half axles and brakes available separately



ITEM #	PRODUCT	SEAT WIDTH CM (IN)	VERSION
530500400	DuoMotion 24 XL	56 (22)	Self Propel
530500500	DuoMotion XL	56 (22)	Transit

Specifications are available on page 68

Item #	OVERALL DIMENSIONS (CM/IN)				SEAT (CM/IN)					WEIGHT (KG/LBS)	
	Width	Length	Height	Height (Front/Rear)	Seat Width	Seat Depth	Seat Height	Opening Width	Opening Length	Product Weight	Weight Capacity
TRS 130 WHEELED COMMODO											
530100000	56 (22)	95 (37)	97 (38)		44 (17)	42,5 (17)	52 (20)			14,4 (31.7)	130 (287)
TS 130 COMMODO											
530100100	58 (23) / 62 (24)	50(20) /55(22)			43 (17)	44 (17)	50-64 (20-25)			8,5 (18.7)	130 (287)
DELUXE DROP ARM COMMODO											
11135-1	71 (28)	47 (19)			58 (23)	47 (19)	48-58 (19-23)			10,0 (22)	300 (661)
BARIATRIC FOLDING COMMODO											
11117N-2	81 (32)		39-56 (15-22)		35 (14)	42 (17)				8,6 (19.0)	295 (650)
TICCO 2G RAISED TOILET SEAT											
523400000	39 (15)	39 (15)		10 (4) / 14.6 (6)				21,5 (8)	28,5 (11)	1,4 (3.0)	225 (496)
523300000	39 (15)	39 (15)		8.6 (3) / 12 (5)				21,5 (8)	28,5 (11)	1,1 (2.3)	225 (496)
TSE 150 RAISED TOILET SEAT											
5201000	56 (22)	48,5 (19)	6,5(2,5)/10(4)/13,5(5,3)					23 (9)	25 (10)	5,7 (12.6)	150 (331)
5200500	46 (18)	48,5 (19)	6,5(2,5)/10(4)/13,5(5,3)					23 (9)	25 (10)	3,7 (8.2)	150 (331)
ELEVATED TOILET SEAT											
12027RA	43 (17)	42 (17)			43 (17)	42 (17)				2,3 (5.1)	135 (298)
COSBY TOILET FRAME											
157	48-61 (19-24)	46 (18)	43-58 (17-23)							4,5 (9.9)	196 (432)
BELLAVITA 2G BATH LIFT											
All items	70 (28)	83 (33)	70 (28)		35 (14)	42 (17)	6,5-48 (3-19)			9,3 (21)	140 (309)
BENNY BATH BENCH											
550100100	69 (27)	35 (14)			40 (16)-63 (25)					2,1 (4.6)	150 (331)
550100200	74,5 (29)	35 (14)			43 (17)-68 (27)					2,6 (5.7)	150 (331)
BEN SWIVELING BATH SEAT											
560600100	73,5 (29")				44 (17)	40 (16)	40 (16)			6,0 (13.2)	120 (265)
DINO SHOWER STOOL											
550500100	45 (18)	45 (18)	52 (20,4)		32 (13)		34-52 (20,4)			1,6 (3.5)	130 (287)
SKY BATH STEP											
RTL12068	50 (20)	40 (16)	10 (4)								140 (309)
DKS 130 FOLDING SHOWER CHAIR											
550700100	48 (19)				41 (16)	41 (16)	46,5-54 (18-21)			4,1 (9.0)	130 (287)
DS 130 SHOWER CHAIR											
550200200	51 (20)	52 (20)			41 (16)	43 (17)	45-53 (18-21)			4,5 (9.9)	130 (287)
DELUXE SHOWER CHAIR											
12009-2	81,3 (32)	51 (20)			68,6 (27)	40,6 (16)	40.6-53.5			6,3 (13.9)	270 (595)
FOXTON STATIONARY SHOWER CHAIR											
140	50 (20)	51 (20)	50 (20)		44 (17)	37 (15)	50 (20)			6,0 (13.2)	165 (364)
UPPINGHAM SHOWER COMMODO											
172	64 (25)	52 (20)			42 (17)	42 (17)	50 (20)			15 (33.1)	160 (353)
ASTON SHOWER COMMODO											
CWC003	56 (22)	90 (35)			46 (18)	43 (17)	51-56 (20-22)			13,4	159 (351)
CWC003HD	66 (26)	90 (35)			56 (22)	43 (17)	51-56 (20-22)			19,4 (42.7)	260 (573)
DUOMOTION SHOWER WHEELCHAIR											
530500100	69 (27)	115 (45)			44 (17)	41 (16)	50-57,5			20 (44.1)	130 (287)
530500200	54,5 (21)	98 (39)			44 (17)	41 (16)	50-57,5			14 (30.9)	130 (287)
DUOMOTION XL SHOWER WHEELCHAIR											
530500400	85 (33)	123 (48)			56 (22)	45 (18)	50-57,5			22,8	200 (441)
530500500	68 (27)	102 (40)			56 (22)	45 (18)	50-57,5			16,8	200 (441)





Kids

Kids Mobility Chairs

Kids Walking Aids

Kids Beds

Kids Bathing Systems



Viper Pediatric Wheelchair

Featuring a lightweight and maneuverable design, the Viper Wheelchair helps to improve activity and mobility. The wheelchair can be folded easily with one hand for storage and transport.

- Available in 12" and 14" seat widths
- Elevating leg rests can be removed
- Padded armrests can pivot backwards

Retractable armrests



Adjustable footrests



Parking brake



ITEM #	SIZE	SEAT WIDTH CM (IN)
L412DDA-ELR	Small	30.5 (12)
L414DDA-ELR	Medium	35.5 (14)

Specifications are available on page 81

Trotter Mobility Chair

The Trotter pediatric stroller offers parents or caregivers a lightweight, foldable alternative to heavier wheelchair-style mobility bases.

- Equipped with 8" front casters and 10" rear wheels with wheel locks
- Adjustable seat depth and swing-away footrests
- Five point harness provides additional support for safety and positioning
- Crash-tested frame for additional vehicle safety

ITEM #	SEAT WIDTH CM (IN)
TR1200	30 (12)
TR1400	35 (14)
TR1600	41 (16)
TR1800	46 (18)

Specifications are available on page 81



Trotter Accessories



Headrest Extension

Provides additional back height and safety during transport.

ITEM# TR8021



Canopy

Height-adjustable and removeable canopy provides sun and rain protection.

ITEM# TR8026



Padded Headrest Wing

Adjustable headrest wing for lateral head positioning and support.

ITEM# TR8022



Lateral Support and Scoli Strap

Provides lateral torso support.

ITEM# TR8027



Utility Bag

Zippered bag for storing personal items.

ITEM# TR8023



Foot and Ankle Positioner

Safely positions and supports feet on footplates.

ITEM# TR8028



Upper Extremity Support Tray

Easy-to-clean and perfect for feeding and interactive play.

ITEM# TR8024



Bus/Transit Tie Downs

Tie down brackets for use in vehicle transport. To comply with safety testing, it should be used with TR8021 & TR8028

ITEM# TR8029



Full Torso Vest

Adjustable vest helps to maintain proper posture and alignment.

ITEM# TR8025

Moxie GT Gait Trainer

The versatile Moxie GT Gait Trainer provides children and adults with stability while they improve their walking ability. This gait trainer can be used in either a posterior or anterior position.

- Height-adjustable in 1" increments and available in three sizes
- Accessories install and adjust without tools



ITEM# GT2000-2GG



ITEM# GT1000-2GB



ITEM# GT3000-2GGY

ITEM #	SIZE	COLOR
GT1000-2GB	Small	Blue
GT1000-2GG	Small	Green
GT2000-2GB	Medium	Blue
GT2000-2GG	Medium	Green
GT3000-2GGY	Large	Grey

Specifications are available on page 81



Moxie Accessories



Moxie Hand Loop

The tool-free hand loops provide stability and can be used as a pair or single.

ITEM# TK1010



Moxie Ankle Prompts

The ankle prompts prevent scissoring of legs and are ideal for users with spasticity, poor motor control, ataxia and excessive adduction.

ITEM# GT1060S / GT1060M / GT1060L



Moxie Guide Bar

Mounted on the center front of the top bar of any size GT, it can be rotated 180° to either pull from the front or push from behind.

ITEM# TK1020



Moxie Hip Positioner

Height adjustable seat harness positions hip and pelvis and allows user to shift weight from one leg to the other

ITEM# TK1070S / TK1070L



Moxie Forearm Platform

Height, depth and angle adjustable for proper arm support.

ITEM# TK1035S / TK1035L



Moxie Trunk Support

Secures the torso at desired angle. Height, depth, width and angle adjustable; Padded for cushioned comfort

ITEM# TK1080S / TK1080M / TK1080L



Moxie Handgrip

Height adjustable and can be mounted in a variety of positions.

ITEM# TK1045



Moxie Thigh Prompt

The height-adjustable thigh prompt limits abduction/adduction and moves with the natural gait pattern.

ITEM# TK1090S / TK1090L

Nimbo Posterior Walker

Lightweight and height-adjustable, the Nimbo Posterior Walker provides users with support for better ambulation. The 5" front and rear rubber wheels can grip any surface without sliding backward.

- Smooth-texture handgrip for added comfort
- Folds easily for transport and storage
- Equipped with override bracket and front caster swivel lock



ITEM# KA2200-2GCR



ITEM# KA1200-2GKB



ITEM# KA4200-2GWP



ITEM# KA5200-2GEB



Nimbo Accessories



Nimbo Forearm Platform

Height, depth and angle adjustable, the forearm platform provides proper arm support.

ITEM# KA1035FPS SIZE Small
ITEM# KA1035FPL SIZE Large



Dial-A-Speed Tabs

Provides variable resistance on one or both rear wheels. They allow user to control gait trainer speed and aligns asymmetrical gait patterns.

ITEM# KA1400 (Pair)



Nimbo Pelvic Stabilizer

The width-adjustable pelvic stabilizer includes a pelvic belt and provides stability to the pelvis for better ambulation.

ITEM# KA1050



Bumper Guards

ITEM# KA8000 (Pair)



Seat Harness

The adjustable harness positions the pelvis and moves with the natural gait pattern of the user.

ITEM# CE1070S SIZE Small
ITEM# CE1070L SIZE Large

NIMBO POSTERIOR WALKER		
ITEM #	SIZE	COLOR
KA1200-2GCR	X-Small	Castle Red
KA1200-2GKB	X-Small	Knight Blue
KA1200-2GWP	X-Small	Wizard Purple
KA2200-2GCR	Small	Castle Red
KA2200-2GKB	Small	Knight Blue
KA2200-2GWP	Small	Wizard Purple

Specifications are available on page 81

NIMBO POSTERIOR WALKER		
ITEM #	SIZE	COLOR
KA3200-2GCR	Medium	Castle Red
KA3200-2GKB	Medium	Knight Blue
KA3200-2GWP	Medium	Wizard Purple
KA4200-2GCR	Large	Castle Red
KA4200-2GKB	Large	Knight Blue
KA4200-2GWP	Large	Wizard Purple
KA5200-2GEB	X-Large	Emperor Black

Specifications are available on page 81



Inspiration Cot

Designed expressly to meet a comprehensive range of pediatric requirements, the Inspiration Cot brings together a range of innovative functions with bright and attractive colors, which combine to make the cot particularly attractive for the pediatric environment.

- Variable height, electrically operated, two-way tilting mattress platform
- Removable head and foot ends allow 360° access when the side rails are down
- Unimpeded access with uniquely designed drop-down telescopic side rails
- Dual sided, linked, tracking and braking system for added safety
- Two removable transfusion poles can be fitted
- Includes specially designed visco elastic mattress for maximum comfort
- Four corner buffers and a shelf for linen
- 9003 model features manual CPR and hinge system on head and foot ends to feed cables/leads and tubes into cot

ITEM# 9000 / 9003

Inspiration Cot Accessories



Traction System

The traction systems are specifically designed for a range of clinical requirements: ORTHOPEDICS, REHABILITATION, NEUROLOGY, EMERGENCY, PEDIATRIC

ITEM# TRACT/PAED



Lullaby Mattress

A mattress developed for use with the Inspiration cot. Manufactured from two foam layers and uses visco elastic foam for the top layer.

ITEM# MAT/LULL/VE

ITEM #	PRODUCT	SPECIFICATION
9000	Inspiration Cot	
9003	Inspiration Cot	with hinged access and CPR
TRACT/PAED	Inspiration Cot Traction System	
MAT/LULL/VE	Lullaby Mattress	



ITEM# OT2000
SIZE Medium

Otter Pediatric Bathing System

Designed with an angle-adjustable seat and backrest to accommodate individual positioning needs, the Otter Pediatric bathing chair provides support and safety for children while bathing.

- Adjustable, slip-resistant legs
- Removable, machine-washable fabric
- Leg straps to control abduction and adduction
- Height- and width-adjustable lateral supports for head and trunk positioning
- Available in three sizes

Adjustable backrest



Adjustable head and neck support



Easy to use



ITEM# OT1000
SIZE Small



ITEM# OT3000
SIZE Large

ITEM #	PRODUCT	SIZE
OT1000	Otter Bathing System	Small
OT2000	Otter Bathing System	Medium
OT3000	Otter Bathing System	Large

Otter Shower Stand

Standing at a comfortable height for the caregiver, the wheeled base of the Otter Shower Stand allows for easy access to roll-in showers.

- Quick-release locks allow bath chair to be easily attached or removed from the frame



Easy to set up and take down

Wheels individually lockable



Otter Tub Stand

With a lightweight, easy-care plastic frame and slip-resistant legs, the Otter Tub Stand adds height to reduce caregiver strain while bathing a child.

- Quick-release locks allow bath chair to be easily attached or removed from frame

2-position height adjustment



ITEM #	PRODUCT	
OT8020	Otter Shower Stand	
OT8010	Otter Tub Stand	

Item #	OVERALL DIMENSIONS CM (IN)				SEAT CM (IN)			ADJUSTMENT	WEIGHT KG (LBS)	
	Width	Length	Height	Folded	Seat Width	Seat Depth	Seat Height	Seat Backrest Adjustment	Product Weight	Weight Capacity
OTTER BATHING SYSTEM										
OT1000			up to 90 (35)	90 x 42,5 x 15 (35 x 17 x 6)	41 (16)	30 (12)	14-26 (6-10)	0-90	4,0 (8.8)	30 (66)
OT2000			80-130 (31-51)	105 x 42,5 x 15 (41 x 17 x 6)	41 (16)	35 (14)	14-26 (6-10)	0-90	4,5 (9.9)	60 (132)
OT3000			110-170 (43-67)	122 x 42,5 x 15 (48 x 17 x 6)	41 (16)	40 (16)	14-26 (6-10)	0-90	5,0 (11.0)	80 (176)
OT8020	68 (27)	90 (35)	51 (20)						4,1 (9.0)	85 (187)
OT8010	42 (17)	78-99 (31-39)	14-27 (6-11)						7,0 (15.4)	85 (187)

Item #	OVERALL DIMENSIONS CM (IN)					SEAT CM (IN)					WEIGHT (KG/LBS)		UOM
	Width	Width Between Arms	Length	Height	Handle Height	Seat Width	Seat Depth	Seat Height	Opening Width	Opening Length	Product Weight	Weight Capacity	
VIPER PEDIATRIC WHEELCHAIR													
L412DDA-ELR													
L414DDA-ELR													
TROTTER MOBILITY CHAIR													
TR1200						30 (12)	30-35	56 (22)			8.6 (19)	34 (75)	
TR1400						35 (14)	33-38	56 (22)			9.5 (21)	45 (99)	
TR1600						41 (16)	35-41	56 (22)			11 (25)	77 (170)	
TR1800						46 (18)	38-43	58 (23)			13.5 (30)	113 (249)	
MOXIE GT GAIT TRAINER													
GT1000-2GB	58 (23)	32 (12.5)	84 (33)		38-47 (15-18.5)						10 (22)	34 (75)	1/cs
GT1000-2GG	58 (23)	32 (12.5)	84 (33)		38-47 (15-18.5)						10 (22)	34 (75)	1/cs
GT2000-2GB	67 (26.5)	42 (16.5)	84 (33)		56-64 (22-25.5)						12 (26)	68 (150)	1/cs
GT2000-2GG	67 (26.5)	42 (16.5)	84 (33)		56-64 (22-25.5)						12 (26)	68 (150)	1/cs
GT3000-2GGY	71 (28)	46 (18)	97 (38)		74-99 (29-39)						15 (32)	91 (200)	1/cs
NIMBO PEDIATRIC WALKING AID													
KA1200-2GCR	61 (24)	34 (13.5)	51 (20)		39-50 (15.5-19.5)						4 (9)	34 (75)	
KA1200-2GKB	61 (24)	34 (13.5)	51 (20)		39-50 (15.5-19.5)						4 (9)	34 (75)	
KA1200-2GWP	61 (24)	34 (13.5)	51 (20)		39-50 (15.5-19.5)						4 (9)	34 (75)	
KA2200-2GCR	61 (24)	34 (13.5)	69 (27)		48-64 (19-25)						5 (10)	39 (85)	
KA2200-2GKB	61 (24)	34 (13.5)	69 (27)		48-64 (19-25)						5 (10)	39 (85)	
KA2200-2GWP	61 (24)	34 (13.5)	69 (27)		48-64 (19-25)						5 (10)	39 (85)	
KA3200-2GCR	66 (26)	38 (15)	79 (31)		58-77 (23-30.5)						6 (13)	68 (150)	
KA3200-2GKB	66 (26)	38 (15)	79 (31)		58-77 (23-30.5)						6 (13)	68 (150)	
KA3200-2GWP	66 (26)	38 (15)	79 (31)		58-77 (23-30.5)						6 (13)	68 (150)	
KA4200-2GCR	67 (26.5)	42 (16.5)	86 (34)		71-91 (28-36)						6 (14)	86 (190)	
KA4200-2GKB	67 (26.5)	42 (16.5)	86 (34)		71-91 (28-36)						6 (14)	86 (190)	
KA4200-2GWP	67 (26.5)	42 (16.5)	86 (34)		71-91 (28-36)						6 (14)	86 (190)	
KA5200-2GEB	74 (29)	48 (19)	94 (37)		91-104 (36-41)						7 (16)	113 (250)	
KA1200S-2GCR	61 (24)	34 (13.5)	51 (20)		39-50 (15.5-19.5)						7 (16)	34 (75)	
KA1200S-2GKB	61 (24)	34 (13.5)	51 (20)		39-50 (15.5-19.5)						7 (16)	34 (75)	
KA1200S-2GWP	61 (24)	34 (13.5)	51 (20)		39-50 (15.5-19.5)						7 (16)	34 (75)	
KA2200S-2GCR	61 (24)	34 (13.5)	69 (27)		48-64 (19-25)						8 (18)	39 (85)	
KA2200S-2GKB	61 (24)	34 (13.5)	69 (27)		48-64 (19-25)						8 (18)	39 (85)	
KA2200S-2GWP	61 (24)	34 (13.5)	69 (27)		48-64 (19-25)						8 (18)	39 (85)	
KA3200S-2GCR	66 (26)	38 (15)	79 (31)		58-77 (23-30.5)						9 (20)	68 (150)	
KA3200S-2GKB	66 (26)	38 (15)	79 (31)		58-77 (23-30.5)						9 (20)	68 (150)	
KA3200S-2GWP	66 (26)	38 (15)	79 (31)		58-77 (23-30.5)						9 (20)	68 (150)	
KA4200S-2GCR	67 (26.5)	42 (16.5)	86 (34)		71-91 (28-36)						19 (21)	86 (190)	
KA4200S-2GKB	67 (26.5)	42 (16.5)	86 (34)		71-91 (28-36)						19 (21)	86 (190)	
KA4200S-2GWP	67 (26.5)	42 (16.5)	86 (34)		71-91 (28-36)						19 (21)	86 (190)	
KA5200S-2GEB	74 (29)	48 (19)	94 (37)		91-104 (36-41)						17 (37)	113 (250)	



▶ **Respiratory**

Oxygen Therapy
Suction Therapy
Aerosol Therapy
CPAP Therapy
Ventilation Therapy
Diagnostics



Compact 525

Delivering up to 5 lpm with high oxygen concentration across all flow rates, this compact concentrator helps to increase patient comfort and lower operating costs.

- Improved sound level reduction and power consumption
- OSD® sensor to continuously monitor oxygen levels
- Visual and audible alarms for low oxygen levels, power failure, pressure drop and service required

ITEM #	PRODUCT
525KS	Compact 525 ultra quiet, International
525KS- LT	Compact 525, International
525PS	Compact 525 ultra quiet, MEA
525DS-Q	Compact 525 ultra quiet, US

Specifications are available on page 100



Compact 1025

Weighing only 19 kg (42 lbs), the 10 Litre Oxygen Concentrator is one of the smallest, most compact high-volume concentrators on the market, perfect for patients who require more than 5 lpm of oxygen. Equipped with exclusive DeVilbiss Oxygen Sensing Device (OSD®) for monitoring oxygen produced, this concentrator helps to ensure user safety and product reliability.

- Capable of delivering up to 96% of O2 purity from 2 to 10 LPM
- Accessible patient controls & clear visible LEDs that can be seen from a distance
- Transfill Caddy (sold separately) with 18 x 14 inch shelf connects easily to 10L Concentrator



ITEM #	PRODUCT
1025KS	Compact 1025

Specifications are available on page 100



iFill Home Oxygen Filling Station

This easy-to-use oxygen filling station extracts and pressurizes room air for quick and convenient cylinder filling. The unlimited, refillable supply of oxygen allows patients to continue filling cylinders as often as required without putting their oxygen therapy on hold.

- Accommodates a range of cylinder sizes for different uses
- Filling station can be placed anywhere in the home
- Comes with unit, mains lead, two cylinders (depending on package ordered) and a cylinder bag



ITEM #	PRODUCT
535I	iFill (Device only)
535I-PD	iFill with 2x PulseDose Cylinders
535I-CF	iFill with 2x Continuous Flow Cylinders
535I-ML6-CF	Continuous Flow Cylinders- ML6
535I-C-CF	Continuous Flow Cylinders- C
535I-D-CF	Continuous Flow Cylinders- D
PD1000A-I-ML6	PulseDose Cylinders- ML6
PD1000A-I-C	PulseDose Cylinders- C
PD1000A-I-D	PulseDose Cylinders- D

Specifications are available on page 100



iGo Portable Oxygen Concentrator

The small and lightweight design of this robust concentrator makes it user-friendly and capable of delivering up to 3 lpm of constant-flow oxygen and 6 lpm (equivalent) in PulseDose mode.

- Supplied with rechargeable lithium-ion battery, AC/DC adapters and mains lead
- Can be operated from standard 12 V car charging socket
- OSD® sensor to continuously monitor oxygen levels



ITEM #	PRODUCT
306DS	iGo (Device only)
306DS-C	iGo with Trolley and Accessory Bag

Specifications are available on page 100



Bonsai Pneumatic Oxygen Conservser

Maintenance-free and lightweight, this oxygen conservser delivers a uniform oxygen pulse in seven dose levels.

- Two continuous flow settings preset at 2 lpm and 4 lpm
- Easy-to-operate single selector knob
- Delivers up to a 6:1 savings ratio from 14-40 breaths per minute
- No batteries required

ITEM #	PRODUCT
OM-812	Connection Type CGA 870 Pin Index
OM-824	Connection Type G 3/4, DIN ISO 477, No.9

Specifications are available on page 87




PD1000 Electronic Oxygen Conservser

This compact oxygen conservser is capable of operating in eight PulseDose mode settings and features adjustable continuous flow flexibility.

- Equipped with protected cannula fitting
- Maximum trigger rate of 40 breaths per minute to avoid hyperventilation
- Acoustic alarm when no breathing is detected
- Built-in pressure indicator and visible battery check
- High-impact polymer case resists damage
- Savings ratio of up to 3:1

ITEM #	PRODUCT
PD1000G	Connection Type G 3/4, DIN ISO 477, No.9
PD1000I	Connection Type CGA 870 Pin Index

Specifications are available on page 87




Evolution Electronic Oxygen Conservser

The lightweight and compact Evolution oxygen conservser delivers a uniform oxygen pulse in seven increments preset at 2 lpm.

- Sensor detects light and heavy exhalation
- Delivers a 5:1 savings ratio from 14-40 breaths per minute
- Minimum 2-year battery life (average 4 hours per day usage) utilizing 2 AA batteries
- Single-push button selects the equivalent lpm dosage settings

ITEM #	PRODUCT
OM-900	Connection Type CGA 870 Pin Index
OM-924	Connection Type G 3/4, DIN ISO 477, No.9
OM-950	Connection Type Quick Connect Fitting
OM-900M	Connection Type CGA 870 Pin Index
OM-924M	Connection Type G 3/4, DIN ISO 477, No.9

Specifications are available on page 87



Item #	OVERALL DIMENSIONS CM (IN)			TECHNICAL DATA			WEIGHT KG (LBS)
	Width	Length	Height	Operating Temperature	Maximum Breathing Rate (pr. minute)	Storing Temperature	Product Weight
BONSAI PNEUMATIC OXYGEN CONSERVER							
OM-812	5.6 (2)	13.5 (5)		-10°C to 40°C	40 breaths	-40°C to 70°C	0.3 (0.6)
OM-824	5.6 (2)	13.5 (5)		-10°C to 40°C	40 breaths	-40°C to 70°C	0.3 (0.6)
ELECTRONIC OXYGEN CONSERVER							
PD1000G		12.4 (4.9)		-10°C to 40°C	40 breaths	-40°C to 70°C	0.6 (1.3)
PD1000I		12.4 (4.9)		-10°C to 40°C	40 breaths	-40°C to 70°C	0.6 (1.3)
EVOLUTION ELECTRONIC OXYGEN CONSERVER							
OM-900	7.9 (3)	15.5 (6)	6.4 (2.5)	-10°C to 40°C	40 breaths	-40°C to 70°C	0.6 (1.3)
OM-924	7.9 (3)	15.5 (6)	6.4 (2.5)	-10°C to 40°C	40 breaths	-40°C to 70°C	0.6 (1.3)
OM-950	7.9 (3)	15.5 (6)	6.4 (2.5)	-10°C to 40°C	40 breaths	-40°C to 70°C	0.6 (1.3)
OM-900M	7.9 (3)	15.5 (6)	6.4 (2.5)	-10°C to 40°C	40 breaths	-40°C to 70°C	0.6 (1.3)
OM-924M	7.9 (3)	15.5 (6)	6.4 (2.5)	-10°C to 40°C	40 breaths	-40°C to 70°C	0.6 (1.3)

VacuAide

The Vacu-Aide Portable Suction Unit is a lightweight, portable aspirator used to clear airways and help control infection.

- Collection bottles include an automatic float shutoff feature to prevent liquid from entering unit
- Can operate on any AC voltage 100 V -240 V
- Comes with carry bag, 100-240 power supply, internal battery, 12 V DC car adapter, mains lead, collection bottle, tubing and filter



ITEM #	PRODUCT
7305P-IS	EU power cord, with battery, 1200ml Autoclavable
7305P-U	UK power cord, with battery, 800ml Disposable
7305P-US	UK power cord, with battery, 1200ml Autoclavable
7305D-G	EU power cord, mains powered only, 800ml Disposable.
7305D-I	EU power cord, mains powered only, 800ml Disposable.
7305D	US power cord, mains powered only, 800ml Disposable

Specifications are available on page 89



VacuAide Compact

This suction unit is compact without compromising suction performance. Featuring a unique float valve that prevents backflow into the unit, and a built-in bacteria filter located in the bottle assembly, this device is optimal for patient care.

- High-performance capabilities with a vacuum range of 50-500 mmHg and a flow range of up to 27 lpm
- Easy-to-read LED bar graph gauge
- ISO 10079-1 conformity

ITEM #	PRODUCT
7310PR-1	VacuAide Compact

Specifications are available on page 89



VacuAide QSU

With a simple on/off operation and large, rotating pressure adjustment gauge, this ultra-quiet unit allows the required suction pressure to be easily set and viewed.

- Features low noise level without compromising suction performance
- Accessories store and transport easily in accompanying carry bag
- ISO 10079-1 conformity

ITEM #	PRODUCT
7314P-NE	EU power cord, With battery, 800ml Disposable
7314P-NE-R	EU power cord, With battery, 1200ml Autoclavable
7314D-NE	EU power cord, Mains powered only, 800ml Disposable
7314D-NE-R	EU power cord, Mains powered only, 1200ml Autoclavable
7314P-U	UK power cord, With battery, 800ml Disposable
7314P-UR	UK power cord, With battery, 1200ml Autoclavable
7314D-U	UK power cord, Mains powered only, 800ml Disposable
7314D-UR	UK power cord, Mains powered only, 1200ml Autoclavable

Specifications are available on page 89

Item #	OVERALL DIMENSIONS CM (IN)			TECHNICAL DATA				WEIGHT KG (LBS)
	Width	Length	Height	Air Flow (l/min)	Vacuum Range (mmHg)	Power	Battery Run Rate (Minutes)	Product Weight
VACUAIDE PORTABLE SUCTION UNIT								
7305P-I	17.8 (7)	20 (8)	22.9 (9)	27	80-550	100 - 240V, 50/60Hz	60	2.9 (6.4)
7305P-IS	17.8 (7)	20 (8)	22.9 (9)	27	80-550	100 - 240V, 50/60Hz	60	2.9 (6.4)
7305P-U	17.8 (7)	20 (8)	22.9 (9)	27	80-550	100 - 240V, 50/60Hz	60	2.9 (6.4)
7305P-US	17.8 (7)	20 (8)	22.9 (9)	27	80-550	100 - 240V, 50/60Hz	60	2.9 (6.4)
VACUAIDE COMPACT PORTABLE SUCTION UNIT								
7310PR-1	18.4 (7)	17 (7)	18.4 (7)	27	50 - 500 ± 10%	100 - 240V, 50/60Hz	45-60	1.53 (3.4)
VACUAIDE QSU PORTABLE SUCTION UNIT								
7314P-NE	20.3 (8)	22 (9)	21.1 (8)	27	50-550 ± 10%	100 - 240 V,50/60 Hz	60	3.0 (6.6)
7314P-NE-R	20.3 (8)	22 (9)	21.1 (8)	27	50-550 ± 10%	100 - 240 V,50/60 Hz	60	3.0 (6.6)
7314D-NE	20.3 (8)	22 (9)	21.1 (8)	27	50-550 ± 10%	100 - 240 V,50/60 Hz	60	3.0 (6.6)
7314D-NE-R	20.3 (8)	22 (9)	21.1 (8)	27	50-550 ± 10%	100 - 240 V,50/60 Hz	60	3.0 (6.6)
7314P-U	20.3 (8)	22 (9)	21.1 (8)	27	50-550 ± 10%	100 - 240 V,50/60 Hz	60	3.0 (6.6)
7314P-UR	20.3 (8)	22 (9)	21.1 (8)	27	50-550 ± 10%	100 - 240 V,50/60 Hz	60	3.0 (6.6)
7314D-U	20.3 (8)	22 (9)	21.1 (8)	27	50-550 ± 10%	100 - 240 V,50/60 Hz	60	3.0 (6.6)
7314D-UR	20.3 (8)	22 (9)	21.1 (8)	27	50-550 ± 10%	100 - 240 V,50/60 Hz	60	3.0 (6.6)



AirForce One Compressor Nebulizer

The AirForce One compact compressor nebulizer offers a cost-effective and reliable form of aerosol therapy. The easy-to-operate and durable unit can be used on adult and pediatric patients.

- Convenient, built-in carry handle
- Comes with necessary accessories and carry bag for immediate use
- 3-year warranty

ITEM #	PRODUCT
MQ5890-EU	AirForce One, EU power cord
MQ5900-UK	AirForce One, UK power cord
<i>Specifications are available on page 100</i>	

Medel Jet Nebulizer



Evo Nebulizer



PulmoAide Compressor Nebulizer

With a carry handle and compact size, the PulmoAide allows for easy transport and effective aerosol therapy.

- Comes with disposable nebulizer
- Easy-to-use design

ITEM #
5650F
<i>Specifications are available on page 100</i>



Reusable Nebulizer



Disposable Nebulizer



UltraNeb Ultrasonic Nebulizer

This high-output, ultrasonic nebulizer is optimal for humidification, sputum removal and aerosol therapy.

- Air flow output and nebulization rates have four adjustable settings and timer function from 0-60 minutes
- Integrated alarm functions alert when chamber is not attached or empty, medication time is complete, or device fails
- Suitable for most sterile water systems
- Compatible with various filling systems

ITEM #	PRODUCT
U3000	Ultrasonic Nebulizer (device only)
U3000H	Ultrasonic Nebulizer with Heated Tube
U3000-DA	Ultrasonic Nebulizer with Desk Model Accessory
U3000H-DA	Ultrasonic Nebulizer with Desk Model Accessory and Heated Tube
U3000-5S	Ultrasonic Nebulizer with 5 Foot Stand
U3000H-5S	Ultrasonic Nebulizer with 5 Foot Stand and Heated tube
<i>Specifications are available on page 100</i>	

Transducer



Front View



Rear View



DeVilbiss Blue StandardPlus

Capable of detecting a wide variety of patient events, the DeVilbiss Blue StandardPlus, a fixed pressure CPAP, is equipped with SmartFlex expiratory relief and Flow Rounding® features.

- SmartCode Rx for remote prescription changes and SmartCode compliance encoding
- Includes SlimFit 15 mm tubing and equipped with MaskFit check
- Easy-to-follow menu, adjustable display brightness and optional clock and alarm clock settings
- Airline transport approved

ITEM #	PAIR W/ HUMIDIFIER (Not Included)
DV63E-NE	DV6HH / DV6HHPD
DV63E-UK	DV6HH / DV6HHPD
DV63E-NOR	DV6HH / DV6HHPD
DV63E-SE	DV6HH / DV6HHPD
<i>Specifications are available on page 101</i>	

DeVilbiss Blue AutoPlus

The DeVilbiss Blue AutoPlus detects, reports and differentiates between hypopneas, obstructive apneas, snoring, expiratory puffs, leaks, flow limitations, RERAs, central apneas and periodic breathing.

- Detects more events and responds with a new level of accuracy
- SmartFlex can be disabled, enabled or used during ramp-time only

ITEM #	PAIR W/ HUMIDIFIER (Not Included)
DV64E-NE	DV6HH / DV6HHPD
DV64E-UK	DV6HH / DV6HHPD
DV64E-NOR	DV6HH / DV6HHPD
DV64E-SE	DV6HH / DV6HHPD
<i>Specifications are available on page 101</i>	



Heated Humidifier:

ITEM #	PRODUCT
DV6HH	DeVilbiss Blue Heated Humidifier
DV6HHPD	DeVilbiss Blue Heated Humidifier with PulseDose
<i>Specifications are available on page 101</i>	

SleepCube Standard

The compact and lightweight SleepCube Standard is a fixed pressure CPAP designed to effectively treat OSA. For added convenience, the heated humidifier easily docks onto the bottom of the unit without requiring separate tubing or power cord.

- Ramp delay settings of 0 to 45 minutes in 5 minute increments
- Automatic altitude compensation up to 2,750 m
- Low-sound pressure level of 26 dBA
- Features auto on/off, visual mask off alert, a large, backlit LCD and easy-to-follow patient menu



ITEM #	PAIR W/ HUMIDIFIER (Not Included)	SMARTCODE TECHNOLOGY (Included)
DV51NE	DV5HH	✓
DV51UK	DV5HH	✓
DV51SE	DV5HH	✓

Specifications are available on page 101

SleepCube AutoPlus

The SleepCube AutoPlus delivers the most effective sleep therapy by using a unique auto-adjusting algorithm that responds to apnea, hypopnea and snoring.

- Equipped with SmartFlex expiratory relief and Flow Rounding comfort features

ITEM #	PAIR W/ HUMIDIFIER (Not Included)	SMARTCODE TECHNOLOGY (Included)
DV54NE-P	DV5HH	✓
DV54UK-P	DV5HH	✓
DV54SE-P	DV5HH	✓

Specifications are available on page 101

SleepCube StandardPlus

With the addition of SmartFlex, an expiratory relief function, the SleepCube StandardPlus helps to lower the average pressure throughout the night for increased user comfort.

- Additional Flow Rounding feature enables a smooth transition from inspiratory to expiratory pressure
- Integrated event detection feature enables the detection of nightly events such as apnea, hypopneas and snoring

ITEM #	PAIR W/ HUMIDIFIER (Not Included)	SMARTCODE TECHNOLOGY (Included)
DV53NE	DV5HH	✓
DV53UK	DV5HH	✓
DV53SE	DV5HH	✓

Specifications are available on page 101

SleepCube AutoBilevel

With an auto-adjusting algorithm that constantly monitors for apneas and events, the SleepCube AutoBilevel adds an extra level of comfort for those who require higher pressures.

- Automatically adjusts to deliver optimal therapy during sleep upon detection of Bilevel pressures

ITEM #	PAIR W/ HUMIDIFIER (Not Included)	SMARTCODE TECHNOLOGY (Included)
DV57NE	DV5HH	✓
DV57UK	DV5HH	✓
DV57SE	DV5HH	✓

Specifications are available on page 101

SleepCube Bilevel S

The SleepCube Bilevel S senses when you are breathing in and out and adjusts the pressure accordingly. This feature is beneficial to users that require higher pressure settings.

- Integrated, real-time event detection feature shows actual pressure, leak rate and patient breathing parameters
- Ten customizable trigger settings for inspiration and expiration

ITEM #	PAIR W/ HUMIDIFIER (Not Included)	SMARTCODE TECHNOLOGY (Included)
DV55NE	DV5HH	✓
DV55UK	DV5HH	✓
DV55SE	DV5HH	✓

Specifications are available on page 101

Heated Humidifier:

ITEM #	PRODUCT
DV5HH	Heated Humidifier

Specifications are available on page 101

SmartLink 3.0 Software Package

The SmartLink software from DeVilbiss Healthcare is a comprehensive therapy management system for monitoring and recording CPAP usage and evaluating the effectiveness of CPAP therapy.



ITEM #	PRODUCT
DV6SS-420	SmartLink 3.0 Software Package

SleepCube Bilevel ST

The SleepCube BiLevel ST gives you the function and performance required to aid in the therapy of respiratory insufficiency, obstructive and central sleep apnoea.

- Integrated, real-time event detection feature shows actual pressure, leak rate and patient breathing parameters
- Ten customizable trigger settings for inspiration and expiration

ITEM #	PAIR W/ HUMIDIFIER (Not Included)	SMARTCODE TECHNOLOGY (Included)
DV56NE	DV5HH	✓
DV56UK	DV5HH	✓
DV65SE	DV5HH	✓

Specifications are available on page 101





J200 Nasal Pillow Mask

The J200 Nasal Pillow Mask offers patients an alternative solution to the conventional CPAP interface. It's small, non-intrusive design minimizes facial contact for optimum comfort. Ideal for patients wearing spectacles during therapy, the mask also offers the patient a larger field of view.

- Discreet and comfortable by design
- Three pillow sizes
- 360° rotating elbow for preferred sleeping position
- Simple ergonomic one-piece design
- Reduced leaking, eliminating skin and eye irritation
- Comfortable headgear including over-head tube positioning
- User friendly, easy to assemble/disassemble for daily cleaning

ITEM# J200



J300 Nasal Mask

Because of its simple design, the lightweight J300 Nasal Mask is perfectly suited for those new to the CPAP therapy. At the same time, experienced CPAP users looking for a solid and reliable nasal mask will value the advantages of the J300. To avoid potential pressure points, the J300 was developed without a front forehead pad.

- Four-component mask enables easy handling and cleaning
- Soft silicon cushion for less leakage
- Contactless separator to the forehead
- Comfortable headband can be changed and adjusted easily

ITEM# J300-M



D100N Nasal Mask

The D100N mask is cost-effective, flexible and easy to maintain. With three sizes of forehead pads and silicone cushions, it is extremely easy to fit the mask and achieve a comfortable seal.

- Fully adjustable headgear
- 360-degree swivel for optimal tubing positioning
- Latex-free



D150F Full Face Mask

Featuring adjustable forehead support and headgear, this mask suits different facial profiles while providing optimal comfort.

- Silicone cushions available in three different sizes: S, M and L
- 360-degree swivel for optimal tubing positioning
- Latex-free



Item #	SIZE	TECHNICAL DATA				WEIGHT G (OZ)
		Storage Temperature Range	Operating Temperature Range	Resistance at 50 l/min	Resistance at 100 l/min	Product Weight (Including Headgear)
J200 NASAL PILLOW MASK						
J200	3 pillows included (size S, M and L)					
J300 NASAL MASK						
J300-M	M					

Item #	SIZE	TECHNICAL DATA				WEIGHT G (OZ)
		Storage Temperature Range	Operating Temperature Range	Resistance at 50 l/min	Resistance at 100 l/min	Product Weight (Including Headgear)
D100N						
D100N-S	S	-40°C to 60°C	-18°C to 50°C	0.47 cmH2O	2.24 cmH2O	88 (3.1)
D100N-M	M	-40°C to 60°C	-18°C to 50°C	0.47 cmH2O	2.24 cmH2O	88 (3.1)
D100N-L	L	-40°C to 60°C	-18°C to 50°C	0.47 cmH2O	2.24 cmH2O	88 (3.1)
D150F						
D150F-S	S	-40°C to 60°C	-18°C to 50°C	0.29 cmH2O	0.86 cmH2O	120 (4.2)
D150F-M	M	-40°C to 60°C	-18°C to 50°C	0.29 cmH2O	0.86 cmH2O	120 (4.2)
D150F-L	L	-40°C to 60°C	-18°C to 50°C	0.29 cmH2O	0.86 cmH2O	120 (4.2)

Cube 30 ATV

A new approach to non-invasive ventilation



The Cube 30 ATV from Drive DeVilbiss International is an innovative, non-invasive volume-supporting ventilator. The Cube 30 ATV incorporates intuitive state-of-the-art functionality, allowing safe and easy application of the device in a homecare or clinical environment. With its unique trigger technology and target minute-volume algorithm, this ventilator optimizes therapy to meet individual patient requirements.

Design

- The device includes a large 4.3" coloured screen display incorporating easy to use tile design menu screens for clear convenient user navigation.



Technology

- High sensitive flow based trigger detection with trigger failure auto-detect function.
- Contemporary volume assured pressure support which can be set based on breath-by-breath target volume or target minute volume.
- Extensive device alarm management, including on screen display alarm protocol.
- Swiss made turbine blower promising high quality performance for extended, reliable operation.

Informative

- Showing clear simple text information the device displays corresponding data settings including user notifications for incorrect parameter settings.
- Three different menu screen colours identifying Patient and Clinical menu modes.



Various Adjustment Options

- Two independent adjustable ventilation settings allowing patients to change therapy mode for daytime and night time use.
- The Mask Check function ensures a correct fitting mask even before therapy starts. On screen values identify leakages in order for the mask to be adjusted and corrected.



Item #	Power Cord	DIMENSIONS	TECHNICAL DATA				WEIGHT KG (LB)
		W x H x D	Storage Temperature Range	Operating Temperature Range	Resistance at 50 l/min	Resistance at 100 l/min	Product Weight
CUBE 30 ATV							
P001S001	EU	16.8x10.8x18.2	-40°C to 60°C	-18°C to 50°C	0.47 cmH2O	2.24 cmH2O	1.6 (3.5)
P001S002	UK	16.8x10.8x18.2	-40°C to 60°C	-18°C to 50°C	0.47 cmH2O	2.24 cmH2O	1.6 (3.5)

SleepDoc Porti 7 Polygraph

The Porti 7 sleep screening device offers 10 standard channels with the option of extending its measurement capabilities up to 24 channels.

The Porti 7 is designed for the diagnosis of Obstructive Sleep Apnoea (OSA) and other sleep related conditions and can be used both in the hospital and at home by the patient.

The device is hard wearing, simple to use and is supplied in a protective carry case.



Nonin WristOX2 Pulse Oximeter

The durable Nonin WristOx2 is ideal for monitoring cardio-ambulatory, remote wireless and overnight studies. This oximeter is compatible with all devices of the DeVilbiss Blue series.

- 48 hours of battery life (non-Bluetooth) and 24 hours with continuous Bluetooth using 2 AAA alkaline batteries
- PureSat signal processing technology
- 1,080 hours of non-volatile memory at a four-second sampling rate with additional high-resolution variable sampling option of one or two seconds
- Convenient data transfer via Bluetooth or USB connection

ITEM #	SPECIFICATION
927701	UK power supply, English Manual
927702	EU power supply, English Manual
927703	EU power supply, French Manual


Specifications are available on page 101



HbO-2000 Fingertip Pulse Oximeter

A non-invasive device designed for oxygen saturation and pulse rate detection, the HbO-2000 is suitable for adult and pediatric patients in the home or clinical setting.

- OLED display includes heart rate bar and plethysmography waveform
- Six different directional screen settings
- Features battery low alert and automatic power off after 8 seconds without detecting a signal
- Low power consumption (2 AAA batteries for 20 hours of operation)



SleepDoc Porti 8 Polysomnograph

The SleepDoc Porti 8 has 10 standard channels and an optional extension of up to 38 channels.

- All SleepDoc Porti systems use the same software, consumables and sensors
- Equipped with a built-in microphone to detect snoring, an integrated magnetic sensor to determine body position and thorax and abdomen belts with pressure sensors
- Capable of programming up to six nights of start/finish times.
- Fully automatic evaluation, diagnosis and report
- Simple operations for both patients and clinical staff
- Fast, direct download via USB connection
- Low-cost consumables reduce running costs
- Internal battery fully charges in 2 hours



ITEM #	SPECIFICATION
928701	UK power supply, English Manual
928702	EU power supply, English Manual
928703	EU power supply, French Manual

Specifications are available on page 101



HbO-3000 Spot Checking Fingertip Pulse Oximeter

Featuring an OLED screen display, SpO2, pulse rate detection and alarms, the stylish, compact HbO-3000 Spot Checking Finger Tip Pulse Oximeter is a must have for patient screening.

- Audible alarms and perfusion index indicator
- Six directional screen display modes and five brightness levels
- Optional plethysmography waveform display
- Includes lanyard and batteries
- Automatic power off

ITEM #	PRODUCT
3150	Nonin WristOx2 Pulse Oximeter
HBO-2000	HbO-2000 Fingertip Pulse Oximeter
HBO-3000	HbO-3000 Spot Checking Fingertip Pulse Oximeter

Specifications are available on page 101

Oxygen Therapy

Item #	OVERALL DIMENSIONS CM (IN)			TECHNICAL DATA				CYLINDER TECH SPECS				WEIGHT KG (LBS)
	Width	Length	Height	Flowrate	Oxygen Concentration	Average Filling	Power	Volume (Litres)	Weight (kg/lbs)	Filling Time	Usage Time	Product Weight
COMPACT 525												
525KS	34.2 (13)	30 (12)	62.2 (24)	0.5-5	93% ± 3%		220-230 V, 50/60 Hz					16.3 (35.9)
525DS-Q	34.2 (13)	30 (12)	62.2 (24)	0.5-5	93% ± 3%		115 V, 60 Hz					16.3 (35.9)
525PS	34.2 (13)	30 (12)	62.2 (24)	0.5-5	93% ± 3%		220-230 V, 60 Hz					16.3 (35.9)
525KS-LT	34.2 (13)	30 (12)	62.2 (24)	0.5-5	93% ± 3%		220-230 V, 50/60 Hz					16.3 (35.9)
COMPACT 1025												
1025KS												
IFILL HOME OXYGEN FILLING STATION												
535I	57 (22)	31 (12)	72.5 (28)		93% ± 3%	138+/- 13.8 bar (2000 +/- 200 psig)	187-254V, 50/60 Hz					30.0 (66.1)
535I-PD	57 (22)	31 (12)	72.5 (28)		93% ± 3%	138+/- 13.8 bar (2000 +/- 200 psig)	187-254V, 50/60 Hz					30.0 (66.1)
535I-CF	57 (22)	31 (12)	72.5 (28)		93% ± 3%	138+/- 13.8 bar (2000 +/- 200 psig)	187-254V, 50/60 Hz					30.0 (66.1)
535I-ML6-CF								1.2	2.18 (5)	90	1.4 h	
535I-C-CF								1.8	2.54 (6)	130	2.0 h	
535I-D-CF								2.9	3.27 (7)	215	3.5 h	
PD1000A-I-ML6								1.2	2.18 (3)	90	4.3 h	
PD1000A-I-C								1.8	2.54 (6)	130	6.1 h	
PD1000A-I-D								2.9	3.27 (7)	215	10.5 h	
IGO PORTABLE OXYGEN CONCENTRATOR												
306DS	27.9 (11)	20 (8)	38.1 (15)	1-3	91% ± 3%							8.6 (19.0), 7.0 (19.0) (w/o battery)
306DS-C	27.9 (11)	20 (8)	38.1 (15)	1-3	91% ± 3%							

Aerosol Therapy

Item #	OVERALL DIMENSIONS CM (IN)			TECHNICAL DATA							WEIGHT KG (LBS)
	Width	Length	Height	Air Flow (l/min)	Particle Size (MMAD)	Supplied Nebulizer (MMAD)	Aerosol Output	Ultrasonic Frequency	Nebulization Rate	Power	Product Weight
AIRFORCE ONE COMPRESSOR NEBULIZER											
MQ5890-EU	14 (6)	19 (7)	9 (4)	5.1	< 5 µm		0.45 ml			230V, 50Hz	1.0 (2.2)
MQ5900-UK	14 (6)	19 (7)	9 (4)	5.1	< 5 µm		0.45 ml			230V, 50Hz	1.0 (2.2)
PULMOAIDE COMPRESSOR NEBULIZER											
5650I	17 (6.5)	26 (10.1)	27 (10.5)	12		< 5 µm			> 0.15	230V, 50Hz	3.2 (7.1)
ULTRANEB ULTRASONIC NEBULIZER											
U3000	20.5 (8)	32 (12)	19 (7)	20	< 4 µm			1.68 MHz ± 5%		230V, 50Hz	3.5 (7.7)
U3000H	20.5 (8)	32 (12)	19 (7)	20	< 4 µm			1.68 MHz ± 5%		230V, 50Hz	3.5 (7.7)
U3000-DA	20.5 (8)	32 (12)	19 (7)	20	< 4 µm			1.68 MHz ± 5%		230V, 50Hz	3.5 (7.7)
U3000H-DA	20.5 (8)	32 (12)	19 (7)	20	< 4 µm			1.68 MHz ± 5%		230V, 50Hz	3.5 (7.7)
U3000-5S	20.5 (8)	32 (12)	19 (7)	20	< 4 µm			1.68 MHz ± 5%		230V, 50Hz	3.5 (7.7)
U3000H-5S	20.5 (8)	32 (12)	19 (7)	20	< 4 µm			1.68 MHz ± 5%		230V, 50Hz	3.5 (7.7)

CPAP Therapy

Item #	OVERALL DIMENSIONS CM (IN)			TECHNICAL DATA						WEIGHT KG (LBS)
	Width	Length	Height	Pressure Range (cmH2O)	Sound Pressure Level	Operating Atmospheric Conditions	Maximum Power Consumption	Average (w/o Humidifier)	Electrical Requirements	Product Weight
DEVILBISS BLUE STANDARDPLUS										
DV63E-NE	15.5 (6.1)	9.4 (3.7)	15 (5.9)	3-20			65 Watts	25 Watts / 10 Watts	100-240 V, 50/60Hz	1.75 (3.8) / 0.9 (1.93)
DV63E-UK	15.5 (6.1)	9.4 (3.7)	15 (5.9)	3-20			65 Watts	25 Watts / 10 Watts	100-240 V, 50/60Hz	1.75 (3.8) / 0.9 (1.93)
DV63E-NOR	15.5 (6.1)	9.4 (3.7)	15 (5.9)	3-20			65 Watts	25 Watts / 10 Watts	100-240 V, 50/60Hz	1.75 (3.8) / 0.9 (1.93)
DV63E-SE	15.5 (6.1)	9.4 (3.7)	15 (5.9)	3-20			65 Watts	25 Watts / 10 Watts	100-240 V, 50/60Hz	1.75 (3.8) / 0.9 (1.93)
DEVILBISS BLUE AUTOPLUS										
DV64E-NE	15.5 (6.1)	9.4 (3.7)	15 (5.9)	3-20			65 Watts	25 Watts / 10 Watts	100-240 V, 50/60Hz	1.75 (3.8) / 0.9 (1.93)
DV64E-UK	15.5 (6.1)	9.4 (3.7)	15 (5.9)	3-20			65 Watts	25 Watts / 10 Watts	100-240 V, 50/60Hz	1.75 (3.8) / 0.9 (1.93)
DV64E-NOR	15.5 (6.1)	9.4 (3.7)	15 (5.9)	3-20			65 Watts	25 Watts / 10 Watts	100-240 V, 50/60Hz	1.75 (3.8) / 0.9 (1.93)
DV64E-SE	15.5 (6.1)	9.4 (3.7)	15 (5.9)	3-20			65 Watts	25 Watts / 10 Watts	100-240 V, 50/60Hz	1.75 (3.8) / 0.9 (1.93)
DEVILBISS BLUE HEATED HUMIDIFIER										
DV6 HHPD										
DV6HH										
SLEEP-CUBE® STANDARD										
DV51NE	16.5 (6.5)	17.5 (7)	10.7 (4.2)	3-25	26 dbA	Sea level to 2600 m (8,500 ft)	65 Watts	28 Watts / 12 Watts	100 - 240 V, 50/60 Hz	1.22 (2.7)
DV51UK	16.5 (6.5)	17.5 (7)	10.7 (4.2)	3-25	26 dbA	Sea level to 2600 m (8,500 ft)	65 Watts	28 Watts / 12 Watts	100 - 240 V, 50/60 Hz	1.22 (2.7)
DV51SE	16.5 (6.5)	17.5 (7)	10.7 (4.2)	3-25	26 dbA	Sea level to 2600 m (8,500 ft)	65 Watts	28 Watts / 12 Watts	100 - 240 V, 50/60 Hz	1.22 (2.7)
SLEEP-CUBE STANDARD PLUS										
DV53NE	16.5 (6.5)	17.5 (7)	10.7 (4.2)	3-25	26 dbA	Sea level to 2600 m (8,500 ft)	65 Watts	28 Watts / 12 Watts	100 - 240 V, 50/60 Hz	1.22 (2.7)
DV53UK	16.5 (6.5)	17.5 (7)	10.7 (4.2)	3-25	26 dbA	Sea level to 2600 m (8,500 ft)	65 Watts	28 Watts / 12 Watts	100 - 240 V, 50/60 Hz	1.22 (2.7)
DV53SE	16.5 (6.5)	17.5 (7)	10.7 (4.2)	3-25	26 dbA	Sea level to 2600 m (8,500 ft)	65 Watts	28 Watts / 12 Watts	100 - 240 V, 50/60 Hz	1.22 (2.7)

Diagnostics

Item #	OVERALL DIMENSIONS (CM/IN)			TECHNICAL DATA					WEIGHT (G/OZ)
	Width	Length	Height	Storage Temperature Range	Operating Temperature Range	SpO2 Measurement Range	Accuracy	Pulse Rate Measurement Range	Product Weight
SLEEPDOC PORTI 7									
927701	6.3 (2)	14 (6)	3.1 (1)						160 (5.6)
927702	6.3 (2)	14 (6)	3.1 (1)						160 (5.6)
927703	6.3 (2)	14 (6)	3.1 (1)						160 (5.6)
SLEEPDOC PORTI 8									
928701	7.5 (3)	17 (7)	3.5 (1)						260 (9.2)
927802	7.5 (3)	17 (7)	3.5 (1)						260 (9.2)
928703	7.5 (3)	17 (7)	3.5 (1)						260 (9.2)
NONIN WRISTOX2									
3150	5 (2)	2 (1)	7 (3)	-40°C to 70°C	-5°C to 40°C	0% to 100%	70% to 100% ± 3%	18 to 321 BPM	
HBO-2000									
HBO-2000	3 (1)	5 (2)	3 (1)	-40°C to 70°C	-5°C to 40°C	35% to 100%	70% to 99% ± 3%	30 to 250 BPM	
HBO-3000									
HBO-3000									

GENERAL TERMS AND CONDITIONS OF DEVILBISS HEALTHCARE GMBH

1. Scope of validity

Sales, deliveries, orders, offers and other services by DeVilbiss Healthcare GmbH (hereinafter referred to as "DeVilbiss") are provided exclusively to business people as defined by Section 14 of the BGB (Bürgerliches Gesetzbuch, German Civil Law) in line with the following General Terms and Conditions that the purchaser recognises by issuing an order or accepting a delivery. They also apply to all future transactions with the purchaser. The applicability of differing or additional terms and conditions issued by the purchaser is excluded even if DeVilbiss does not expressly contradict them.

2. Agreeing the contract

- 2.1. DeVilbiss offers are non-binding. A contract only comes into existence when DeVilbiss accepts the purchaser's order.
- 2.2. The purchaser may issue orders via the internet, fax or by telephone.
- 2.3. After receiving the purchaser's order over the internet DeVilbiss will send a non-binding confirmation of receipt to the purchaser.
- 2.4. DeVilbiss will accept the purchaser's order by sending an order confirmation to the purchaser by fax or email. The contract depends exclusively on the content of the order confirmation and the present General Terms and Conditions. Oral subsidiary agreements or commitments must be confirmed in writing by DeVilbiss for them to be effective.
- 2.5. DeVilbiss reserves all rights to the sales documents (especially graphics, drawings, and weight and measurement statements) and samples. They may not be made accessible to third parties and must be returned to DeVilbiss without delay on request. Such documents may only be passed on to third parties after the prior written agreement of DeVilbiss.
- 2.6. The DeVilbiss sales team cannot sign contracts of any kind nor make binding statements on the item to be delivered or other conditions.

3. Delivery periods and deadlines

- 3.1. Delivery dates and periods are only binding if they have been confirmed by DeVilbiss in writing and the purchaser has provided or made available to DeVilbiss all the information and documents required to execute the delivery and has paid any advance payments as agreed. Agreed periods start on the date of order confirmation. The periods are extended as appropriate for additional or extension orders issued later.
- 3.2. Unpredictable, unforeseeable events as well as those beyond the influence or responsibility of DeVilbiss, such as force majeure, war, natural disasters, release DeVilbiss from the duties of on-time delivery or service provision for their duration. The agreed periods are extended for the duration of the disruption; the purchaser will be informed of the occurrence of the disruption in an appropriate manner. If the end of the disruption can not be predicted or it lasts longer than two months, each party is entitled to withdraw from the contract.
- 3.3. If the purchaser delays acceptance or infringes its other cooperation obligations, notwithstanding its other rights, DeVilbiss is entitled to store the goods at the purchaser's expense and risk as appropriate.
- 3.4. DeVilbiss may make partial deliveries if there is a good reason to do so.

4. Shipment, transfer of risk and insurance policies

- 4.1. If no provision is made by the purchaser, shipment is made in the usual packaging by an appropriate delivery method.
- 4.2. Shipment takes place ex works, which means that DeVilbiss' liability ends when the consignment is handed over to the transport company or, if the purchaser is responsible for transport, when the consignment is provided at the DeVilbiss warehouse from which shipment takes place and the announcement of readiness for shipment. In this case, even if the purchaser is not responsible for any damages to the consignment or its loss, the purchaser is still obligated to pay DeVilbiss the full amount of the purchase price.
- 4.3. Insurance policies are only provided at the request and the expense of the purchaser.
- 4.4. For all orders that fall below an order value of €50 excl. VAT, DeVilbiss reserves the right to add a packing and handling surcharge of €15 excl. VAT.

5. Prices and payment terms

- 5.1. If the ordering parties have not agreed on a particular price, the price is set as per the DeVilbiss price list applicable at the time the contract is concluded.
- 5.2. All DeVilbiss prices are understood as ex works excluding the relevant, applicable value-added tax, public fees and charges, any customs charges as well as packaging and shipment costs, which are charged separately.
- 5.3. The parties agree that if nothing different has been negotiated; all payment obligations arising from this contract are in euros.
- 5.4. If nothing further has been agreed the invoice sum is to be paid by direct debit (cash in advance) immediately after delivery against 3% cash discount. If payment on invoice has been agreed each invoice is due for payment without deduction immediately after receipt. Payments by the purchaser are considered to have been made when DeVilbiss can dispose of the amount.
- 5.5. If the purchaser has delayed payment DeVilbiss is entitled to demand default interest as permitted by law. DeVilbiss may demand restitution of any other damage due to delay without limitation.
- 5.6. The purchaser is only entitled to offset amounts if the counter-claim is undisputed or legally binding.
- 5.7. Purchaser is only permitted to execute a retention right to the extent to which the counter-claim is based on the same contract and is undisputed or legally binding.
- 5.8. If DeVilbiss becomes aware of the risk that the purchaser can not fulfil his duties after the contract has been concluded, DeVilbiss is entitled to only make outstanding deliveries or provide services only for payment in advance or if security is provided; if this is not done even after a period set by DeVilbiss, DeVilbiss may withdraw from the contract in full or in part notwithstanding other rights.

6. Warranty, inspection requirement and return

- 6.1. Data given in catalogues, price lists and other handouts supplied to the purchaser by DeVilbiss and in product descriptions should in no way be taken as a guarantee that the article shall have a particular characteristic. Any such guarantee must be agreed explicitly in writing.
- 6.2. The purchaser's warranty rights require that the purchaser shall check the article following delivery and that DeVilbiss shall be notified immediately in writing of any defect by no later than 2 weeks from the delivery date. Latent defects must be notified to DeVilbiss in writing immediately upon discovering the defect. All warranty rights are not transferable to subsequent owners.
- 6.3. DeVilbiss reserves the right to inspect and test the rejected article in any claim. The purchaser will give DeVilbiss the necessary time and opportunity to do so. DeVilbiss may also request after arrangement that the purchaser returns the item complained about to DeVilbiss at the expense of DeVilbiss. If a complaint by the purchaser proves to be unjustified due to deliberate action or negligence, DeVilbiss is entitled to receive reimbursement for all costs incurred in this regard.
- 6.4. DeVilbiss shall be entitled, if it so chooses, to remedy any warranty defects by repairing or replacing the defective part or the entire article at no extra cost to the purchaser.

GENERAL TERMS AND CONDITIONS OF DEVILBISS HEALTHCARE GMBH

- 6.5. Any parts replaced by DeVilbiss will become the property of DeVilbiss.
- 6.6. DeVilbiss does not accept any responsibility for damages arising from inappropriate or incorrect use, stocking, storage, transportation, assembly, operation, and lack of maintenance or inappropriate handling by the purchaser, use of inappropriate fixtures or natural wear and tear, except where those damages are caused by DeVilbiss.
- 6.7. If a warranty defect is not remedied by rectification or replacement or if such is unacceptable to the purchaser, or if DeVilbiss has rejected doing so by virtue of Article 439, Section 3 of the BGB, the purchaser may choose to cancel the agreement in respect of the defective article, request a corresponding reduction in price or request simple compensation, compensation in lieu of specific performance or reimbursement of its futile expenditures.
- 6.8. The limitation period for warranty claim is 12 months from the transfer of risk. If the item to be delivered is sold properly by the purchaser or their direct or indirect contractual partners to an end-user, the provisions of Section 479 of the BGB for any recourse claims remain in force for a longer limitation period.
- 6.9. The purchaser shall not be entitled to any other claim except where allowed under Section 7 or under Articles 478 and 479 of the BGB.
- 6.10. Finished products or parts thereof cannot be returned if there are no grounds for so doing. If in individual cases, especially for therapeutic reasons, modified requirements or deaths DeVilbiss takes a different decision, DeVilbiss will pay 80% of the delivery price for original packed goods that were delivered within the last three months. Deliveries that were made more than 3 months ago are excluded from exchange. A copy of the delivery note or invoice is proof of the delivery date. Freight costs are to be assumed by the purchaser.

7. Damages and liability restriction

- 7.1. If nothing further has been agreed in Number 7.2 DeVilbiss is liable without limitation for damages as per the statutory regulations.
- 7.2. As an exception DeVilbiss is not liable for slight negligence of insignificant contractual obligations and only to a limited extent for damages that can be predicted as being typical of the contract for slight negligence of key contractual obligations.

8. Reservation of title


- 8.1. The products delivered remain the property of DeVilbiss until all claims by DeVilbiss from the business relationship with the purchaser have been paid in full.
- 8.2. In case of current account title retention serves to guarantee the balance claim assignable to DeVilbiss.
- 8.3. The purchaser is only permitted to sell the retained products within the ordinary course of business. The purchaser is not permitted to pledge the retained products, transfer them as security or otherwise use them in a manner that could put DeVilbiss ownership at risk. The purchaser already assigns the claim from ongoing sales to DeVilbiss; DeVilbiss already accepts this assignment. The purchaser is revocably authorised to collect the claims assigned to DeVilbiss in its own name in trust for DeVilbiss. DeVilbiss can revoke this authorisation and the right to sell the goods on if the purchaser is in default with key obligations such as payment.
- 8.4. The purchaser will provide to DeVilbiss at any time all of the desired information on the retained products or claims that are hereby assigned to DeVilbiss. The purchaser must notify DeVilbiss immediately of accesses to or claims by third parties on the retained products and provide the necessary documents. The purchaser will also notify the third parties of the DeVilbiss retention of title. The costs for defending such accesses and claims will be paid by the purchaser.
- 8.5. The purchaser is obliged to treat the retained products carefully for the duration of the retention of title.
- 8.6. If the realisable value of the securities exceeds the total protected claims of DeVilbiss by more than 10% the purchaser is entitled to request release to this extent.
- 8.7. If the purchaser is delayed in key duties, such as payment to DeVilbiss, DeVilbiss can notwithstanding other rights remove the retained products and for the purposes of fulfilling any claims against the purchaser make use of them in another way. In this case the purchaser will grant DeVilbiss or those appointed by DeVilbiss immediate access to the retained products and hand them over. If DeVilbiss demands the return as a result of this provision this does not mean withdrawal from the contract.
- 8.8. For deliveries in other legal systems where the retention of title provision stated above does not have the same security effect as in Germany, the purchaser will do everything in its power to give DeVilbiss appropriate security rights without delay. The purchaser will cooperate on all measures, such as registration, publication etc. that are necessary and appropriate for the effectiveness and ability to implement such security rights.
- 8.9. At the request of DeVilbiss the purchaser is obliged to insure the retained products appropriately, to provide appropriate confirmation of the insurance and to assign claims from the insurance policy to DeVilbiss.

9. Product liability / medical products

- 9.1. The purchaser will not change the delivered products nor their equipment or packaging; in particular he will not modify or remove the existing warnings about risks relating to the improper use of the products. If the purchaser infringes the provision stated above he will indemnify DeVilbiss internally for product liability claims by third parties if he is responsible for the errors that result in liability.
- 9.2. Obligation to cooperate in corrective measures according to Medical Products Safety Scheme Directive (MPSV).
- 9.3. The purchaser is obligated to carry out, forward and report back about carrying out corrective measures which had been communicated by DeVilbiss.

10. General provisions

- 10.1. Modifications and additions to the contract and/or these general terms and conditions, as well as subsidiary agreements, must be in written form. This also applies to any changes to this written form clause.
- 10.2. If a provision of the contract and/or these general terms and conditions is ineffective in part or in full this does not affect the effectiveness of the other provisions. In this case the parties agree to be obliged to agree in the place of the ineffective provision on a provision that comes as close as possible to the economic purpose of this ineffective provision.
- 10.3. The exclusive place of jurisdiction for all disputes arising from the present contractual relationship is Mannheim. DeVilbiss is however entitled to issue proceedings in any other statutory court of jurisdiction.
- 10.4. These terms and conditions are subject to the laws of the Federal Republic of Germany with the exclusion of the United Nations Convention on Contracts for the International Sale of Goods (CISG).
- 10.5. In order to conduct business processes your information will be stored and processed by DeVilbiss. DeVilbiss will not transmit your information to any other third party. As long as you don't raise an explicit or in any other way to DeVilbiss discernible objection, DeVilbiss will use your information for marketing purposes.



UK Office

DeVilbiss Healthcare Ltd
Unit 3 Bloomfield Park
Bloomfield Road
Tipton, West Midlands
DY4 9AP
Phone: + 44 (0) 121 521 3140
Fax: + 44 (0) 121 521 3141
eMail: uk@drivedevilbiss-int.com
Open: Monday to Friday 09.00 – 17.00

German Office

DeVilbiss Healthcare GmbH
Kamenzer Str. 3
68309 Mannheim
Germany
Phone: + 49 (0) 621 178 98 - 0
Fax: + 49 (0) 621 178 98 -111
eMail: germany@drivedevilbiss-int.com
Open: Monday to Thursday 08.30 – 17.00
Friday 08.30 – 15.00

Indian Office

DeVilbiss Healthcare pvt Ltd.
Unit no.18
1st Floor Shree Madhu Estate
P.B. Marg, Worli
Mumbai 400013
India
Phone: +91 22 24960000
eMail: india@drivedevilbiss-int.com
Open: Monday to Friday
09.30 –18.00

Australian Office

Drive DeVilbiss pyt Ltd.
Unit 8, 15 Carrington Road
Castle Hill,
New South Wales 2154
Australia
Phone: + 61 2 9899 3144
Fax: + 61 2 9899 3244
eMail: australia@drivedevilbiss-int.com.au
Open: Monday to Friday 09.00 – 17.00



DeVilbiss Healthcare GmbH • Kamenzer Str. 3 • 68309 Mannheim • Germany
Phone: +49 (621) 178 98 0 • Fax: +49 (621) 178 98 111 • kontakt@devilbisshc.com • www.drivedevilbiss-int.com

November 2017

Drive DeVilbiss Healthcare - All content, images, trademarks and trade names in association with this catalogue are owned by Drive DeVilbiss Healthcare. The information in this catalogue has been provided by Drive DeVilbiss Healthcare for general information purposes only. While we endeavour to keep the information up to date and correct; we make no representations or warranties of any kind about the accuracy, reliability and suitability of the information. No part of this catalogue may be reproduced, adapted or copied in any form without prior written permission from Drive DeVilbiss Healthcare. Permission is only deemed valid if approval is in writing.

drive

DeVilbiss
HEALTHCARE